IPS e.max® ZirCAD labside

[en] Instructions for Use

- Zirconium oxide disc (ZrO₂)
- Zirconium oxide block (ZrO₂)
- Colouring Liquid

[de] Gebrauchsinformation

- Zirkoniumoxid-Scheibe (ZrO₂);
- Zirkoniumoxid-Block (ZrO₂)
- Färbeliquid

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Patients with permanent teeth

Intended users / Special training

- Dentists (chairside fabrication of restorations; clinical workflow)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)

No special training required.

Use

For dental use only.

Description

- IPS e.max® ZirCAD is a yttrium-stabilized zirconium oxide for fixed all-ceramic dental restorations.
- IPS e.max ZirCAD Colouring Liquids are used for shading and characterizing IPS e.max ZirCAD restorations.

Technical data

	IPS e.max ZirCAD Prime Esthetic		IPS e.max ZirCAD Prime		IPS e.max ZirCAD MT Multi	
Shade range*	BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4		BL1, BL2, BL3, BL4, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4		BL1, BL3, A1, A2, A3, A3.5, B1, B2, C1, C2, C3, D2, D3	
Disc thicknesses	14 mm, 16 mm, 20 mm		14 mm, 16 mm, 20 mm, 25 mm		16 mm, 20 mm	
Block sizes	–		–		C17, B45	
	Specifications	Typical average value	Specifications	Typical average value	Specifications	Typical average value
Flexural strength [MPa]	≥ 700	850	≥ 900	1200	≥ 700	850
Fracture toughness	3.6 MPa • m ^{1/2} (Dentin)		>5.0 MPa • m ^{1/2} (Dentin)		3.6 MPa • m ^{1/2} (Dentin)	
CTE (25 – 500°C) [10 ⁻⁶ K]	10.4 ± 0.5		10.5 ± 0.5		10.4 ± 0.5	
Chemical solubility [µg/cm²]	< 100		< 100		< 100	
MDR classification	Class IIa		Class IIa		Class IIa	
Type / Class ISO 6872:2015	Type II / Class 4		Type II / Class 5		Type II / Class 4	
	IPS e.max ZirCAD MT		IPS e.max ZirCAD LT		IPS e.max ZirCAD MO	
Shade range*	BL, A1, A2, A3, B1, B2, C2, D2		0, 1, 2, 3, 4, sun, sun chroma BL, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4		0, 1, 2, 3, 4	
Disc thicknesses	14 mm, 18 mm		10 mm, 12 mm, 14 mm, 16 mm, 18 mm, 20 mm, 25 mm		10 mm, 14 mm, 18 mm, 20 mm, 25 mm	
Block sizes	–		C17, B45		C13, C15, C15L, B40, B40L, B55, B65, B65L-17, B85L-22	
	Specifications	Typical average value	Specifications	Typical average value	Specifications	Typical average value
Flexural strength [MPa]	≥ 700	850	≥ 900	1200	≥ 900	1150
Fracture toughness	3.6 MPa • m ^{1/2}		5.1 MPa • m ^{1/2}		5.1 MPa • m ^{1/2}	
CTE (25 – 500°C) [10 ⁻⁶ K]	10.4 ± 0.5		10.5 ± 0.5		10.5 ± 0.5	
Chemical solubility [µg/cm²]	< 100		< 100		< 100	
MDR classification	Class IIa		Class IIa		Class IIa	
Type / Class ISO 6872:2015	Type II / Class 4		Type II / Class 5		Type II / Class 5	

Indications

- **IPS e.max ZirCAD MO, LT, Prime (incl. IPS e.max ZirCAD LT Colouring Liquids)**
Missing tooth structure in anterior and posterior teeth, partial edentulism in the anterior and posterior region (max. 2 pontics)
- **IPS e.max ZirCAD MT, MT Multi, Prime Esthetic (incl. IPS e.max ZirCAD MT Colouring Liquids)**
Missing tooth structure in anterior and posterior teeth, partial edentulism in the anterior and posterior region (max. 1 pontic, max. 3-unit bridges)

Types of restorations

Anterior and posterior restorations on prepared teeth and on approved implant abutment systems

	Full-contour crowns	Full-contour 3-unit bridges	Full-contour 4- and multi-unit bridges with max. 2 pontics	Crown copings	3- and multi-unit bridge frameworks with max. 2 pontics
IPS e.max ZirCAD Prime	✓	✓	✓	✓	✓
IPS e.max ZirCAD Prime Esthetic	✓	✓			
IPS e.max ZirCAD MT Multi	✓	✓			
IPS e.max ZirCAD MT	✓	✓			
IPS e.max ZirCAD LT	✓	✓	✓	✓	✓
IPS e.max ZirCAD MO				✓	✓

Contraindications

- Patients with substantially reduced residual dentition
- The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients

Limitations of use

- Any other uses that are not included in the range of applications.
- Temporary insertion.
- Untreated bruxism (a splint is indicated after incorporation).
-  Do not reuse. The final restoration must not be reused.

Processing restrictions

Failure to observe the following restrictions may compromise the success achieved:

- Failure to observe the necessary minimum wall thicknesses and connector dimensions
- Milling the discs in a non-compatible CAD/CAM system
- Sintering in a non-compatible sintering furnace

System requirements

A milling system that is suitable for the processing of 98.5-mm discs must be used.

Side effects

There are no known side effects to date.

Interactions

There are no known interactions to date.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Dental ceramics	IPS e.max ZirCAD Prime Esthetic	IPS e.max ZirCAD Prime	IPS e.max ZirCAD MT Multi	IPS e.max ZirCAD MT	IPS e.max ZirCAD LT	IPS e.max ZirCAD MO
Zirconium oxide (ZrO ₂)	86.0 – 93.5%	87.0 – 95.5%	86.0 – 93.5%	86.0 – 93.5%	88.0 – 95.5%	88.0 – 95.5%
Yttrium oxide (Y ₂ O ₃)	>6.5 – ≤ 8.3%	> 4.5 – ≤ 7.0%	> 6.5 – ≤ 8.0%	> 6.5 – ≤ 8.0%	> 4.5 – ≤ 6.0%	> 4.5 – ≤ 6.0%
Hafnium oxide (HfO ₂)	≤ 5.0%	≤ 5.0%	≤ 5.0%	≤ 5.0%	≤ 5.0%	≤ 5.0%
Aluminium oxide (Al ₂ O ₃)	≤ 1.0%	≤ 1.0%	≤ 1.0%	≤ 1.0%	≤ 1.0%	≤ 1.0%
other oxides	≤ 1.0%	≤ 1.5%	≤ 1.0%	≤ 1.0%	≤ 1.0%	≤ 1.0%

2 Application

Procedure:

- CAD design
- CAM nesting
- Milling
- Finishing
 - Optional: Brush infiltration + drying
- Sintering
- Finishing
 - Optional: Veneering technique
- Staining and glazing

Application notes

I. Minimum thicknesses and connector dimensions

IPS e.max ZircAD Prime Esthetic	Types of restorations	Anterior region		Posterior region	
		Minimum layer thickness in mm	Connector dimensions** in mm ²	Minimum layer thickness in mm	Connector dimensions** in mm ²
	Crowns	0.8	–	1.0	–
	3-unit bridges	1.0	12*	1.0	16

* Height: 4 mm, width: 3 mm

** The specified minimum connector cross-section must be positioned in the dentin area of the disc

IPS e.max ZircAD Prime	Types of restorations	Anterior region		Posterior region	
		Minimum layer thickness in mm	Connector dimensions** in mm ²	Minimum layer thickness in mm	Connector dimensions** in mm ²
	Fully or partially anatomical restorations				
	Crowns	0.8	–	1.0	–
	3-unit bridges	1.0	9	1.0	9
	4- and multi-unit bridges with 2 pontics***	1.0	9	1.0	12*
	Cantilever bridges with one pontic	1.0	12*	1.0	12*
	Frameworks; prerequisite: complete positioning in the dentin area				
	Crowns	0.4	–	0.6	–
	3-unit bridges	0.6	9	0.6	9
	4- and multi-unit bridges with 2 pontics***	0.6	9	1.0	12*
	Cantilever bridges with one pontic	1.0	12*	1.0	12*

* Height: 4 mm, width: 3 mm

** The specified minimum connector cross-section must be positioned in the dentin area of the disc

*** In Canada, bridge indications are limited to 6 units with a maximum of 2 connected pontics

IPS e.max ZircAD MO / LT	Types of restorations	Anterior region		Posterior region	
		Minimum layer thickness in mm	Connector dimensions in mm ²	Minimum layer thickness in mm	Connector dimensions in mm ²
	Crowns	0.4	–	0.6	–
	3-unit bridges	0.6	7	0.6	9
	4- and multi-unit bridges with 2 pontics**	0.6	9	0.7	12*
	Cantilever bridges with one pontic	0.7	12*	0.7	12*

* Height: 4 mm, width: 3 mm

** In Canada, bridge indications are limited to 6 units with a maximum of 2 connected pontics


IPS e.max ZircAD MT / MT Multi	Types of restorations	Anterior region		Posterior region	
		Minimum layer thickness in mm	Connector dimensions** in mm ²	Minimum layer thickness in mm	Connector dimensions** in mm ²
	Crowns	0.8	–	1.0	–
	3-unit bridges	1.0	12*	1.0	16

* Height: 4 mm, width: 3 mm


** IPS e.max ZircAD MT Multi: The specified minimum connector cross-section must be positioned in the dentin area of the disc

II. Positioning guidelines

When using IPS e.max ZirCAD Prime and Prime Esthetic

When using IPS e.max ZirCAD MT Multi

General notes

- The height positioning in the polychromatic disc determines the amount of incisal area in the restoration.
- Given the esthetic characteristics of single-tooth restorations, it is advisable to place the restoration about 1 mm below the top edge of the disc, in order to ensure a clearly visible incisal area.
- For fully or partially anatomical restorations, make sure that the specified minimum connector cross-section is positioned in the dentin area of the disc.
- Crown and bridge frameworks must be placed completely in the dentin area.

III. Nesting

Detailed information can be found in the Instructions for Use of the CAM software.

General notes

- 3 holding bars must be attached to each single-tooth restoration.
- In multi-unit restorations, end units must have 2 holding bars attached (oral and vestibular). Attach holding bars to other units as required.
- The diameter of the holding bars must be at least 2.0 mm.
- The holding bars must be attached at least 1.0 mm above the preparation edge.
- The holding bars should be placed at the anatomical equator of the tooth so that no undercuts are created and the restoration can be easily processed from both sides.
- In the case of multi-unit restorations with a pronounced curvature, a sintering support structure is advisable.
- Provide end units with a holding bar vertical to the sintering support structure.
- Design the sintering support structure in an even thickness (2 – 5 mm).

IV. Milling

- The disc information can be transferred to any compatible CAM software that is able to read the RFID label.
- For manual entry the specific shrinkage factor can be found on the disc. Detailed information can be found in the Instructions for Use of the CAM software.
- The imprinted side corresponds with the incisal / occlusal area.
- When securing the disc in the holder, care must be taken to ensure that the circumferential notch and the disc holder are absolutely clean and that the screws are tightened evenly crosswise.

V. Separating and finishing

- Suitable rotary tools (e.g. fine tungsten carbide burs) are recommended for separating the restorations. It is advisable to notch the holding bars on one side before the restoration is completely separated.
- Suitable rotary tools (e.g. fine tungsten carbide burs) are recommended to smooth out the attachment points of the holding bars. Rough tungsten carbide burs and/or grinding instruments are not suitable, as they may cause chipping, among other things.
- When using the brush infiltration technique, make sure the surfaces of the restorations are not contaminated with any grease or smudgy film since this may negatively influence the shade result.
- After finishing, thoroughly remove any zirconium oxide dust with a soft brush and then blast the restoration with oil-free compressed air.

General notes

- Non-sintered zirconium oxide restorations are susceptible to damage and fractures. This fact has to be kept in mind during the entire working procedure.
- Any adjustments should be carried out in the non-sintered state.
- In the non-sintered state, any contact with unsuitable liquids and liquids not approved for zirconium oxide (e.g. unpurified water and/or lubricant coolant) and/or contact media (e.g. occlusion spray) must be prevented.
- Use only light pressure for finishing.
- Do not, under any circumstances, post-separate bridge constructions with a separating disc. This may lead to predetermined breaking points in the area of the connectors and reduce the strength of the all-ceramic restoration.
- Rubber polishers should not be used as they condense the surface and cause contamination.
- Make sure that the minimum wall thicknesses of the restorations are maintained during finishing.
- Adhering zirconium oxide dust may be sintered to the restoration and lead to fitting inaccuracies.
- The non-sintered restoration must not be sonicated or steam-cleaned.
- The non-sintered restoration must not be blasted.

Brush infiltration

The following aspects should be taken into consideration when applying the brush infiltration technique:

- The IPS e.max ZirCAD MT Colouring Liquids and IPS e.max ZirCAD LT Colouring Liquids must not be mixed or combined. They may only be used on the intended zirconium oxide.
- The colouring liquids have been specially developed for the brush infiltration technique and must not be used for the dip infiltration technique.
- The restoration must be free from dust and grinding residue.
- Restorations must be completely dry before infiltration.

- The colouring liquids must not be contaminated.
- The colouring liquids have to be sealed when not in use.
- If there is cloudiness or precipitation (e.g. sediments), the colouring liquids should no longer be used. Cloudiness or precipitation of the colouring liquids is promoted by contamination.
- Do not decant and/or store the colouring liquids in metal containers.
- Apply IPS e.max ZirCAD Colouring Liquids with a clean, metal-free brush on the restoration.
- Do not use IPS e.max ZirCAD Colouring Liquid Indicator in an unmixed state.
- IPS e.max ZirCAD Colouring Liquid Indicator colours are not stable for longer periods of time.
- Store mixed solutions in a sealed container and use them within 4 hours. After longer storage, the reproducibility of the tooth shade in the sintered state is no longer ensured.
- Infiltrated restorations must be completely dried before sintering.

Warnings

It is advisable to wear protective gloves during the infiltration process. They help prevent skin irritation caused by the colouring liquids and the formation of a grease film on the restoration, which may compromise the infiltration of the colouring liquids.

VIII. Sintering

Sintering is a thermal process during which IPS e.max ZirCAD obtains its final physical and mechanical properties, such as its high strength and translucency.

In general, the following notes should be taken into consideration when sintering:

- Only fully dried restorations may be sintered in order to prevent damage to the furnace and/or the object.
- The use of sintering beads is not recommended.
- Sufficient atmospheric exchange in the sintering sagger must be ensured.
- The restorations must not come into contact with each other during sintering.
- The correct program selection must be observed.
- Sintering temperatures that are too low or too high and/or sintering times that are too short or too long will have a negative effect on the above mentioned final properties.
- Always keep the sintering accessories clean and free of dust so that the sintered restorations are not contaminated.

Information on sintering

Standard program for the sintering of IPS e.max ZirCAD Prime/Prime Esthetic/MT Multi/MT/LT/MO crowns / bridges with up to 14 units and infiltrated, dried restorations.

Standard program up to 14 units	Temperature 1 [°C/°F]	Temperature 2 [°C/°F]	Heating rate [°C/°F/min]	Holding time [min]
Heating phase	20/68	900/1652	10/18	–
Holding phase	900/1652	900/1652	–	30
Heating phase	900/1652	1500/2732	3/6	–
Holding phase	1500/2732	1500/2732	–	120
Cooling phase	1500/2732	900/1652	10/18	–
Cooling phase	900/1652	300/572	8/15	–
– switch off –				

Speed program for the quick sintering of IPS e.max ZirCAD Prime/Prime Esthetic/MT Multi/MT/LT/MO crowns / bridges with up to 3 units.

Speed program up to 3 units	Temperature 1 [°C/°F]	Temperature 2 [°C/°F]	Heating rate [°C/°F/min]	Holding time [min]
Heating phase	20/68	1000/1832	60/108	–
Holding phase	1000/1832	1000/1832	–	10
Heating phase	1000/1832	1530/2786	3/6	–
Holding phase	1530/2786	1530/2786	–	60
Cooling phase	1530/2786	1100/2012	50/90	–
Cooling phase	1100/2012	100/212	60/108	–
– switch off –				

Processing after sintering

The following aspects should be taken into consideration when processing restorations after sintering:

- Processing of sintered restorations should be kept to a minimum.
- The restoration should only be mechanically processed if it is absolutely necessary.
- Use only light pressure and low speed to adjust the restoration.
- Prevent sharp edges when finishing frameworks.
- Bridge connectors must not be post-separated with a separating disc.
- We recommend using a rubber polisher to smooth the basal side of the bridge connectors.
- The material-specific minimum wall thicknesses and connector dimensions must not be undermined during processing.
- Use only immaculate grinding instruments
- Remove zirconium oxide dust sintered to the restoration with suitable grinding instruments. As an alternative, zirconium oxide dust sintered to the restoration may be removed by blasting with Al_2O_3 25 – 70 µm at a pressure of 1 bar (15 psi) or 70 – 110 µm at a pressure of 1.5 bar (22 psi).

- Rinse the restoration under running water or use the steam jet to remove any adhering residue and dry.
- Finish the restoration in accordance with the desired processing technique (staining, cut-back or layering).

Preparing for cementation

Condition the IPS e.max ZirCAD restoration by sandblasting the inner aspects of the crown with Al_2O_3 , 25 – 70 μm , 1 bar (15 psi) or Al_2O_3 , 70 – 110 μm , 1.5 bar (22 psi).

3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).
- The Summary of Safety and Clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed>. Basic UDI-DI: 76152082ACERA006F2

Warnings

- Observe the safety data sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).
- The processing of IPS e.max ZirCAD discs and blocks produces dust which may irritate the skin and eyes and which may result in lung damage. Make sure that the suction equipment at your milling machine and at your workplace works flawlessly. Do not inhale grinding dust during finishing and wear a dust mask (particle class FFP2) as well as protective goggles and gloves.

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Chipping / fracture / decementation of the restorative material may lead to accidental swallowing or inhalation and dental re-treatment.
- Cement excess may lead to irritation of the soft tissue / gingiva. Progressive inflammation may lead to bone resorption or periodontitis.

4 Shelf life and storage

- In the original packaging
- In a dry place
- Do not expose to physical impact or vibrations.

IPS e.max ZirCAD Colouring Liquids:

- Store at 2 – 28 °C / 36 – 82 °F.
- Do not use the product after the indicated expiry date.
- Expiry date: see note on bottles and packages.
- Before use, visually inspect the packaging and the product for damage. In case of any doubt, consult Ivoclar Vivadent AG or your local trade partner.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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