IPS InLine®

[en] Instructions for Use

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- Patients with permanent teeth

Intended users

- Dentists (clinical procedure)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)

Special training

None

Use

For dental use only.

Description

IPS InLine is a leucite-containing veneering ceramic suitable for the fabrication of metal-ceramic restorations at firing temperatures above $900\text{ }^{\circ}\text{C}$ ($1652\text{ }^{\circ}\text{F}$).

IPS InLine has been classified as a Type 1, Class 1 dental ceramic.

IPS InLine can be used for veneering of alloys in the CTE range of $13.8\text{--}15.2 \times 10^{-6}/\text{K}$ ($25\text{--}500\text{ }^{\circ}\text{C}$) and for fabricating veneers on refractory dies.

Product name	Product description
IPS InLine System Opaquer IPS InLine Gingiva Opaquer IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	Masking of metal frameworks and adjustment of the base shade
IPS InLine Margin	Fabrication of ceramic margins on metal frameworks
IPS InLine Deep Dentin IPS InLine Dentin	Veneering of metal frameworks; fabrication of layered veneers in the anterior region
IPS InLine Incisal IPS InLine Transpa Incisal	Veneering of the incisal and occlusal portions of metal frameworks; fabrication of layered veneers in the anterior region
IPS InLine Transpa	Characterization of the veneering of metal frameworks; Characterization of layered veneers in the anterior region
IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon IPS InLine Opal Effect IPS InLine Cervical Incisal	Characterization of the veneering of metal frameworks; Characterization of layered veneers in the anterior region
IPS InLine One Denticisal	Veneering of metal frameworks
IPS InLine Add-On Margin IPS InLine Add-On IPS InLine System Add-On $690\text{ }^{\circ}\text{C}/1274\text{ }^{\circ}\text{F}$	Adjustment of ceramic margins and restorations on metal frameworks
IPS InLine Gingiva	Veneering of metal frameworks in the gingival area
IPS InLine Intensive Gingiva	Characterization of Gingiva materials

Accessories for the medical device

Product name	Product description
IPS InLine System Opaquer Liquid	For diluting IPS InLine System Opaquer
IPS InLine System Powder Opaquer Liquid	For mixing IPS InLine System Powder Opaquer
IPS InLine Margin Build-Up Liquid	For mixing IPS InLine Margin layering materials
IPS InLine System Build-Up Liquid (L, P)	For mixing IPS InLine layering materials

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Types of restorations:

- Crowns
- 3–14-unit bridges
- Veneers

Contraindications

- The use of this product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

The product must not be used in the following cases:

- Untreated bruxism
- Mixing and processing together with other dental ceramics
- Failure to observe the necessary minimum connector and framework thicknesses
- Exceeding or falling short of the stipulated veneering layer thicknesses
- Failure to observe the layer thickness ratio between the framework and layering ceramics
- Veneering of dental alloys not within the stipulated CTE range
- Veneering of titanium and zirconium oxide frameworks
- Reuse of the final restoration

Side effects

None known to date

Interactions

None known to date

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Product components powders and pastes		Ingredients
Opaquers	IPS InLine System Opaquer IPS InLine Gingiva Opaquer	– Leucite glass-ceramic – Glass – Zirconium oxide – Inorganic pigments – Glycols (butandiol, glycerine) – Polyvinylpyrrolidone – Highly dispersed silicic acid
	IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	– Leucite glass-ceramic – Zirconium oxide – Inorganic pigments
Margin	IPS InLine Margin	– Leucite glass-ceramic – Inorganic pigments
Deep Dentin Dentocisal Dentin Incisal Transpa Incisal	IPS InLine Deep Dentin IPS InLine One Dentocisal * IPS InLine Dentin * IPS InLine Incisal * IPS InLine Transpa Incisal *	– Leucite glass-ceramic – Opal glass – Inorganic pigments
Impulse	IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensive Gingiva IPS InLine Opal Effect	
Add-On	IPS InLine Add-On	– Opalescent leucite glass-ceramic – Glass – Inorganic pigments
	IPS InLine Add-On Margin	– Leucite glass-ceramic – Inorganic pigments
	IPS InLine System Add-On 690 °C/1274 °F	– Leucite glass-ceramic – Glass – Inorganic pigments

*) additionally contains small amounts of identification colour, which burns without residue when the ceramic is fired.

Product components ceramic mixing liquids	Ingredients
IPS InLine System Opaquer Liquid	– Glycols (butandiol, glycerine) – Polyvinylpyrrolidone
IPS InLine System Powder Opaquer Liquid	– Water – Glycol (propandiol) – Metallic salt – Carbonic acid

Product components ceramic mixing liquids	Ingredients
IPS InLine Margin Build-Up Liquid	– Water – Cellulose derivative
IPS InLine System Build-Up Liquid (L, P)	– Water – Polyethylene glycol – Metallic salt

2 Application

NOTICE! Observe the following instructions:

- Wall thickness for metal frameworks: Crowns min. 0.3 mm; abutment crown min. 0.5 mm.
- Ceramic layering thickness: min. 0.8 mm; max. 1.7 mm.

NOTICE! Only use the system components mentioned in chapter "Intended use". Proceed as follows:

- Mix powder materials with mixing liquid before use.
- Extrude and mix paste materials before use. Use mixing liquid for diluting.

NOTICE! Conduct all firing cycles according to the instructions of the manufacturer. Firing parameters can be found in chapter Technical specifications.

2.1 Layering crowns/bridges

Fabricating the framework

1. Fabricate and oxidize the framework with a reduced anatomical shape according to the instructions of the alloy manufacturer.

Applying and firing the opaquer

1. Apply a thin, non-covering layer of paste or powder opaquer on the framework.
2. Conduct first opaquer firing (wash firing).
3. Optional: For reinforced in-depth fluorescence, mix 20% of IPS InLine System Opaquer F with the paste opaquer.
4. Apply a thin, covering layer of paste or powder opaquer on the framework.
5. Apply an even covering layer of Gingiva paste or powder opaquer on the gingiva portions
6. Customize the framework using Intensive paste or powder opaquer.
7. Conduct second opaquer firing.
8. Optional: For reinforced in-depth fluorescence apply IPS InLine System Opaquer F in a thin third opaquer layer.
9. Optional: Conduct third opaquer firing (firing parameters of 1/2. paste opaquer firing)

A One-layer technique

Layering and firing One material

1. Fully build up the restoration using the one-layer material
2. Conduct One firing.
3. Repeat the steps 1 to 2.

B Standard and individual layering technique

In case of ceramic margins: Layering and firing the ceramic margin

1. Isolate and seal the die in the margin areas.
2. Apply Margin material in the cervical area
3. Dry the Margin material.
4. Lift the framework with Margin material off the die.
5. Conduct Margin firing.
6. Check the fit.
7. Remove interfering areas.
8. Repeat the steps 2–7.

Layering and firing Dentin and Incisal material

1. If space is limited: Apply opaque Dentin material on the opaquerized framework.
2. Completely layer the restoration with Dentin material.
3. Cut back the restoration.
4. Add Incisal material to the incisal area.
5. Optional: Customize the restoration using characterization material.
6. Conduct Dentin and Incisal firing.
7. Repeat the steps 4 to 6.

Layering and firing Gingiva material

1. Layer the gingiva area.
2. Optional: Customize the restoration using characterization material.
3. Conduct Gingiva firing.
4. Repeat the steps 1 to 3.

Characterizing, glazing and firing restorations with IPS Ivocolor

NOTICE! Make sure that the adhesive surface of the restoration is completely free of ceramic material.

1. Finalize the anatomical shape and surface structure.
2. Customize the restoration using stains.
3. Conduct stain firing.
4. Apply glaze material using a brush.
5. Conduct glaze firing.

Applying and firing corrections

NOTICE! Corrections can be fired with the glaze firing cycle or after the glaze firing cycle.

1. **Correction with glaze firing:** Apply the Add-On material.
2. Conduct glaze firing.
3. **Correction after glaze firing:** Apply the Add-On 690 °C/1274 °F material.
4. Conduct corrective firing.

2.2 Layering veneers on refractory dies

Duplicating the die

1. Duplicate the die in silicone.
2. Remove the die.
3. Pour a commercial refractory die material into the silicone mould.
4. Degas the refractory die.
5. Soak the die in water.

Applying and firing the wash

1. Mix IPS InLine Add-On with IPS Ivocolor Mixing Liquid allround.
2. Apply the material in a thin layer on the refractory die.
3. Conduct wash firing.

Layering and firing the veneer

1. Layer Dentin and Impulse material on the die.
2. Conduct Dentin/Impulse firing.
3. Build up incisal areas using Incisal materials.
4. Conduct Incisal firing.

Characterizing, glazing and firing the restoration

See chapter Layering crowns/bridges.

Divesting the veneer

1. Remove the die material from the veneer using a polishing jet medium (max. 1 bar/15 psi).

Conditioning the veneer for adhesive cementation

1. Blast the inner aspect of the veneer with Al₂O₃ (max. 1 bar/15 psi).
2. Clean the veneer with the steam jet.
3. Dry the veneer with oil-free air.
4. Etch the inner aspect of the veneer with IPS Ceramic Etching Gel for 60 s.
5. Repeat the steps 2 to 3.

2.3 Seating the restoration

Seating the restoration

NOTICE! Seating is exclusively done by the dentist.

- Seat crowns and bridges using conventional or adhesive methods.
- Seat veneers using adhesive methods.

3 Technical specifications

Physical-chemical properties of product components of IPS InLine according to 1) ISO 6872:2024 and 2) ISO 9693:2019.

Property / Product components	Radioactivity [Bq/g U ²³⁸]		Biaxial flexural strength σ [MPa]		Chemical solubility [µg/cm ²]		Adhesive bond Tb [MPa]	
	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	C ³⁾	B ²⁾
IPS InLine System Opaquer IPS InLine Gingiva Opaquer	≤ 1,0	<< 1,0	> 50	170	< 100	28	> 25	40
IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	≤ 1,0	< 0,03	> 50	127	< 100	28	> 25	43
IPS InLine Margin	≤ 1,0	<< 1,0	> 50	98	< 100	23	N/A	N/A
IPS InLine Deep Dentin IPS InLine One Dentsisal IPS InLine Dentin IPS InLine Incisal IPS InLine Transpa Incisal	≤ 1,0	<< 1,0	> 50	89	< 100	12	> 25	40
IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensive Gingiva	≤ 1,0	<< 1,0	> 50	89	< 100	12	> 25	40
IPS InLine Opal Effect	≤ 1,0	<< 1,0	> 50	90	< 100	10	N/A	N/A
IPS InLine Add-On	≤ 1,0	<< 1,0	> 50	92	< 100	10	N/A	N/A

Property / Product components	Radioactivity [Bq/g U ²³⁸]		Biaxial flexural strength σ [MPa]		Chemical solubility [μ g/cm ²]		Adhesive bond Tb [MPa]	
	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	A ¹⁾	B ²⁾	C ³⁾	B ²⁾
IPS InLine Add-On Margin	≤ 1,0	<< 1,0	> 50	94	< 100	12	N/A	N/A
IPS InLine System Add-On 690 °C/1274 °F)	≤ 1,0	<< 1,0	> 50	108	< 100	20	N/A	N/A

1) A: Acceptance criterion according to ISO 6872:2024

2) B: Typical measurements (mean values)

3) C: Acceptance criterion according to ISO 9693:2019

Thermal properties of product components of IPS InLine according to 1) ISO 6872:2024 and 2) ISO 9693:2019.

Properties / Product components	Coefficient of linear thermal expansion (CTE)			Glass transition temperature t_g [°C]	Thermal cycling test
	CTE 2x ¹⁾ [10 ⁻⁶ /K]	CTE 4x ²⁾ [10 ⁻⁶ /K]	CTE $\bar{\sigma}$ ³⁾ [10 ⁻⁶ /K]		
IPS InLine System Opaquer IPS InLine Gingiva Opaquer	13,5 ± 0,5	13,7 ± 0,5	13,6 ± 0,5	605 ± 20	Applicable
IPS InLine System Powder Opaquer IPS InLine System Intensive Powder Opaquer IPS InLine Gingiva Powder Opaquer	13,2 ± 0,5	13,4 ± 0,5	13,3 ± 0,5	605 ± 20	Applicable
IPS InLine Margin	13,2 ± 0,5	13,5 ± 0,5	13,4 ± 0,5	600 ± 20	Applicable
IPS InLine Deep Dentin IPS InLine One Denticisal	12,4 ± 0,5	13,2 ± 0,5	12,8 ± 0,5	580 ± 20	N/A
IPS InLine Dentin IPS InLine Incisal IPS InLine Transpa Incisal	12,4 ± 0,5	13,2 ± 0,5	12,8 ± 0,5	580 ± 20	Applicable
IPS InLine Cervical Incisal IPS InLine Cervical Dentin IPS InLine Occlusal Dentin IPS InLine Mamelon IPS InLine Transpa IPS InLine Gingiva IPS InLine Intensive Gingiva	12,4 ± 0,5	13,2 ± 0,5	12,8 ± 0,5	580 ± 20	N/A
IPS InLine Opal Effect	12,6 ± 0,5	13,5 ± 0,5	13,1 ± 0,5	595 ± 20	N/A
IPS InLine Add-On	12,3 ± 0,5	12,8 ± 0,5	12,6 ± 0,5	455 ± 20	N/A
IPS InLine Add-On Margin	13,3 ± 0,5	13,7 ± 0,5	13,5 ± 0,5	585 ± 20	N/A
IPS InLine System Add-On 690 °C/1274 °F *)	12,6 ± 0,5	13,3 ± 0,5	13,0 ± 0,5	440 ± 20	N/A

1) CTE after two firing cycles, measurement range 25–500 °C/77–932 °F

2) CTE after four firing cycles, measurement range 25–500 °C/77–932 °F

3) CTE mean value (after two/four firing cycles), measurement range 25–500 °C/77–932 °F

*) The CTE values refer to a temperature range of 25–400 °C/77–752 °F, as the glass transition temperature is below 500 °C/932 °F.

Firing table

NOTICE! These firing parameters are guidance values and apply to Programat® furnaces from Ivoclar Vivadent AG. Deviations (approx. ± 10 °C/ 50 °F) may occur:

- depending on the furnace generation
- with ceramic furnaces from other manufacturers
- due to regional differences in the power supply
- when operating several electrical devices in the same circuit

IPS InLine	Firing temperature T [°C]	Stand-by temperature B [°C]	Closing time S [min]	Heating rate t [°C/min]	Holding time H [min]	Vacuum on V1 [°C]	Vacuum off V2 [°C]
1./2. Powder opaquer firing	960/1760	403/757	4:00	100/180	2:00	450/842	959/1758
1./2. Paste opaquer firing	930/1706	403/757	6:00	100/180	2:00	450/842	929/1704
1./2. Margin firing	930/1706	403/757	4:00	60/108	1:00	450/842	929/1704
1. Dentin / Incisal / Gingiva / Dentical firing	910/1670	403/757	4:00	60/108	1:00	450/842	909/1668
2. Dentin / Incisal / Gingiva / Dentical firing	900/1652	403/757	4:00	60/108	1:00	450/842	899/1650
Margin Add-On firing	930/1706	403/757	4:00	60/108	1:00	450/842	929/1704
Corrective firing after Dentin / Incisal firing, Add-On	860/1580	403/757	4:00	60/108	1:00	450/842	859/1578
Stain firing using IPS Ivocolor	830/1526	430/757	6:00	60/108	1:00	450/842	829/1524
Glaze firing using IPS Ivocolor	830/1526	403/757	6:00	60/108	1:00	450/842	829/1524
Add-On after glaze firing (690 °C/1274 °F)	690/1274	403/757	4:00	60/108	1:00	450/842	689/1274

IPS InLine Veneer	Firing temperature T [°C]	Stand-by temperature B [°C]	Closing time S [min]	Heating rate t [°C/min]	Holding time H [min]	Vacuum on V1 [°C]	Vacuum off V2 [°C]
Wash firing	830/1526	403/757	4:00	60/108	1:00	450/842	829/1524
Cervical firing	940/1724	403/757	8:00	60/108	1:00	450/842	939/1722
Dentin / Impulse firing	940/1724	403/757	8:00	60/108	1:00	450/842	939/1722
Incisal firing	930/1706	403/757	8:00	60/108	1:00	450/842	929/1704
Stain/glaze firing using IPS Ivocolor materials	830/1526	403/757	8:00	60/108	1:00	450/842	829/1524

4 Safety information

- The product has been developed solely for use in dentistry. Process according to the Instructions for Use.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- Do not inhale ceramic dust during finishing. Use suction equipment, a face mask and protective goggles.
- If applied frequently or for a prolonged period of time, peroxide-based professional bleaching agents (carbamide peroxide; hydrogen peroxide) as well as acidulated phosphate fluorides used for caries prevention can make the surface of restorations rough and matte.
- The product contains hexane. Hexane is highly flammable and detrimental to health. Avoid contact of the material with the skin and eyes. Do not inhale vapours. Keep material away from sources of ignition.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

Supporting documents

Document	To be found at:
Current version of Instructions for Use	www.ivoclar.com/eIFU
Structure of Instructions for Use and Warnings	www.ivoclar.com/eIFU
Explanation of symbols	www.ivoclar.com/eIFU
Safety Data Sheet (SDS)	www.ivoclar.com
Summary of Safety and Clinical Performance – SSCP	https://ec.europa.eu/tools/eudamed
Basic-UDI-DI: 76152082ACERA011ET	

Disposal information

Remaining stocks or removed restorations must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Chipping / fracture / debonding of the restorative material may lead to accidental swallowing or inhalation and dental retreatment

5 Shelf life and storage

- Storage temperature ceramic pastes and mixing liquids: 2–28 °C (36–82 °F)
- Storage temperature for ceramic powders: No special storage conditions are required.
- Date of expiration: See note on bottles, syringes or packaging
- Store the product free of vibration.
- Store the product in a dry place.
- Protect the product from sunlight.

6 Additional information

Keep material out of the reach of children!

The manufacturer accepts no liability for damages resulting from improper use of the product. Furthermore, the user is responsible for determining the product's suitability for the user's own purposes before use.

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