






Garrison



Product Catalogue
2024



MATRIX SYSTEMS

Strata-G™ 4-5

Deep Margin Elevation 22-23 & 25

Reel Matrix™ 24-25

Slick Bands™ Tofflemire-style Matrices 26-27

Fusion™ Anterior 30-31

VariStrip™ Anterior 29

Blue View™ Cervical Matrices 35

RINGS

Strata-G™ 6

3D Fusion™, 3DXR, 3D500, Gold, Original 7

BANDS

Strata-G™ 9

Firm, Fusion, SXR 10

SM, Gold, M, B 11

WEDGES

Strata-G™ 12

Fusion, A+, G-Wedge™, WedgeWands™ 13

PRE-WEDGE SHIELD

FenderWedge® 15

INTERPROXIMAL FINISHERS

FitStrip™ 16-19

POLISHING SYSTEMS

eZr™ 20-21

MATRIX BAND FORCEPS

Composi-Tight® Matrix 34

RING PLACEMENT FORCEPS

Premier 35

Original 35

Universal 35

CONTACT FORMERS

Perform™ 33

COMPOSITE INSTRUMENT

Multi-Function Non-Stick 36

Standard Non-Stick 36

INFORMATION PAGES

Strata-G™ 2-4

Cases 3, 8, 15, 16, 22, 23, 28, 30

CE Information 14

Matrix Selection Guide 38-39

Order Form 40

Ordering Information inside back cover

Guarantee back cover

Garrison® sectional matrix systems are the standard of care.

Tofflemire-style System


- Fails to restore proximal anatomy
- Thin contact at the marginal ridge
- Large food trap (below)
- Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease



Strata-G™ Sectional Matrix System

- + Operator-friendly retaining system
- + Naturally contoured bands
- + Anatomically correct contacts
- + Contacts at the height of contour
- + Rings produce optimal tooth separation for excellent, tight contacts

Now faster and more predictable!



Reliable, tight contacts

Drawn-wire nickel titanium plus PEEK produces consistent separating pressure longer!

Easier, more predictable contours

Marginal ridge enhancements built into the matrix bands AND rings help provide that finishing flair!

Faster, less frustrating Class IIs

The Strata-G™ rings, bands and wedges combine for a 270° Seal™ for ultimate flash control. The rings won't spring off the tooth, and the amazing Wide Preparation ring makes tough wide preps no problem.

Garrison® sectional matrix systems are the standard of care.


Garrison® Sectional Matrix Cases

Great adaptation to the teeth eliminates the need for adjustment with polishers in the proximal areas.



Images courtesy of Dr. Javid Muraslov, Latvia

Wide preps can be handled in one simple step when you have the right (green) ring.



Images courtesy of Dr. Clarence Tam, HBSc, DDS, FIADFE, AAACD, Auckland, New Zealand

Back-to-back requires a shift in technique—we can show you how.*



Images courtesy of Dr. Anand Narvekar, Mumbai, India

Quadrants can be a challenge. Ease your concerns with tips and tricks from the pros.



Images courtesy of Prof. Nicola Scotti, DDS, PhD Torino, Italy

The Award Winning Sectional Matrix System



TESTIMONIALS

David Naegele, DMD
(The Strata-G™ ring)
created a much better seal on the prep for moisture control, and minimal, if any, flash. Strata-G™ bands produced clear marginal ridges as I've never seen before. Finishing time for the restoration is GREATLY reduced.

Nick Romenesko, DDS
The rubber feet of the (Strata-G™) rings help reduce flash, create a better seal around the box and my marginal ridge roll is excellent. Since receiving Strata-G™ it is the only system I use... I would cry if someone took it away.

Steven Berwitz, DMD

I have used various sectional matrix systems from Garrison®, and I like this one the best. There is a significant reduction in flash and the contours of the restorations are ideal. The contacts made by this system are great!

*To see full writeups, cases and reviews go to www.garrisondental.com/learning
Do you have a case you want to share? Email it to us at contact@garrisondental.com.

STRATA-G™


Strong tooth separation yet easier to open!

Nitinol drawn-wire nickel titanium is aligned at the molecular level for superior performance yet requires less effort to open. Garrison's proprietary StrataSet™ treatment maximizes the longevity of this advanced material.

More uses!

The PEEK ultra-durable reinforcement on the back of the ring allows Strata-G™ to retain a higher percentage of initial strength for much longer than non-reinforced rings. Proprietary StrataBond™ silicone/PEEK tips withstand hundreds of applications.


Improved anatomy—Unbeatable grip



Wider indications—wider smiles!

Conserve more of your patient's natural tooth structure with dentistry's premier Wide Preparation ring. The Strata-G™ Wide Preparation ring (green) provides a more conservative alternative for severely broken-down teeth.

Strata-G™ forms a 270° Seal™ for ultimate flash control.



STRATA-G™

faster—more predictable—flash control



SG-KSH-11 kit contains:

- 1 Universal/Short Ring (blue)
- 1 Tall Ring (orange)
- 1 Wide Prep Ring (green)
- 50 Assorted Matrix Bands
- 50 Firm Matrix Bands
- 10 VariStrips™
- 80 Assorted Ultra-Adaptive Wedges
- 5 deep margin elevation bands
- 1 Ring Placement Forceps (pg 35)
- 1 BFA-1 Band Placement Forceps (pg 34)
- 1 TN009 Multi-Function Instrument (pg 36)

Item

SG-KSH-11

SG-KSH-10 (BFA-1 and TN009 not in kit)

SG-KSH-14 (TN009 not in kit)

For refill information on Strata-G™ rings, bands, and wedges, refer to pages 6, 9 and 12.

For more kit configurations and systems, go to garrisondental.com

STRATA-G™

Longer lasting by design.

Advanced materials and design combine to produce our longest lasting, most predictable separator ring ever!



All-New Separator Rings.

- The Strata-G™ Universal/Short ring (blue) can easily grip the distal of the canine, short or malpositioned teeth and even primary teeth. The Universal ring can help you restore the widest range of cases in dentistry.
- The Strata-G™ Tall ring (orange) provides ideal matrix band adaptation and separation on taller teeth. Additionally, it is perfect for stacking over the Universal ring (blue) for MODs or multiple tooth restorations.





Dentistry's Premier Wide Preparation Ring.

- The Strata-G™ Wide Preparation ring (green) allows you to conserve more of your patient's natural tooth structure. It provides a more conservative alternative than typical indirect restorations for severely broken-down teeth. Treat these cases as an extra-wide Class II rather than a complete crown prep with the confidence of knowing that you will achieve excellent contact and contour.

STRATA-G™

Smaller, easier to open and lasts longer

		
SG400 Universal/Short Ring 2/pk	SG500 Tall Ring 2/pk	SG600 Wide Prep Ring 2/pk

Composi-Tight. 3DFusion.

Soft-Face™ silicone tips

		
FX400 Short Ring 2/pk	FX500 Tall Ring 2/pk	FX600 Wide Prep Ring 2/pk

Composi-Tight. 3D XR

The tenacious grip, along with enhanced Soft-Face™ technology and Dynamic Tip Angle™


3DXR Soft-Face™ Extra Retention Rings 2/pk



Composi-Tight. 3D

The Soft-Face™ 3D Ring provides maximum adaptation to the tooth


3D500 Soft-Face™ Rings 2/pk

Composi-TightGold.

2x stronger than our Original G-Rings™

	
AU400 Standard Length Tines 3/pk	AU500 Long Length Tines 3/pk

Composi-Tight.

The original G-Ring™ that started it all—made of polished stainless-steel

	
G100 Standard Length Tines 3/pk	G200 Long Length Tines 3/pk



Case using Firm Matrix Bands

Dr. Anand Narvekar (Mumbai, India) utilized Garrison® Firm Bands for deep margin elevation in this challenging onlay case.




Figure 1. FXH300 large molar band with extension was placed. The wedge was wrapped with multiple layers of teflon tape to sufficiently increase its size to allow it to seal the cervical portion of the matrix band in this extremely wide interproximal space.



Figure 2. Cervical area was built up with multiple layers of composite.



Figure 3. The process was repeated for the adjacent molar. The stiff nature of Firm Bands allows apical pressure to be applied to the band without deformation of the matrix. It can be seated deeper than dead soft matrices for these types of procedures.





Figure 4. One-week post-bonding of onlays.

STRATA-G™

Anatomy top to bottom

Five sizes of super-thin stainless steel in both dead-soft and firm ensure you have the best matrix band to restore every case.



Optimized marginal ridge

Strata-G™ Bands provide shape and guidance for developing the marginal ridge of Class II restorations without blocking your access.

ProfileBooster™ subgingival apron

Achieving a natural emergence profile in challenging cervical areas is easy with Strata-G™ Bands—it's built right in!



- Non-stick and color coded
- Easy placement tabs
- Dead-soft and only 0.0016" thin
- NEW pedodontic size!

Strata-G™ Matrix Bands

Description	Item	Qty.
Assorted Band Kit: 5 sizes	SGB01	400
Assorted Band Kit: 5 sizes	SGB02	200






SG Series

Strata-G™	Strata-G™ matrix bands with optimized marginal ridge				
	Pediatric	Premolars	Premolars and small molars	Molars	Molars with deep cervical preps
thickness 0.0016"					
height	SG050 3.5 mm	SG100 4.5 mm	SG175 5.5 mm	SG200 6.5 mm	SG300 7.5 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
SG050	50	100
SG100	50	100
SG175	50	100
SG200	50	100
SG300	50	100

Band Comparison Chart

FXH Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement.						
thickness 0.0015"						
height	NA	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm

Item	Standard (-M) Qty.	Bulk Qty.
FXH100	50	100
FXH150	50	100
FXH175	50	100
FXH200	50	100
FXH300	50	100

FX Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.						
thickness 0.0016"						
height	NA	FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm	FX300 6.0 mm w/ext 8.7 mm

Item	Standard (-M) Qty.	Bulk Qty.
FX100	50	100
FX150	30	60
FX175	50	100
FX200	50	100
FX300	30	60

SXR Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.						
thickness 0.0016"						
height	NA	SXR100 4.6 mm	SXR150 3.8 mm w/ext 5.4 mm	SXR175 5.4 mm	SXR200 6.4 mm	SXR300 6.4 mm w/ext 8.6 mm

Item	Standard (-M) Qty.	Bulk Qty.
SXR100	50	100
SXR150	30	60
SXR175	50	100
SXR200	50	100
SXR300	30	60

To order standard packs, add a -M after the item #.

Band Comparison Chart

SM Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Slick Bands™ Slick Bands™ Matrices have non-stick coating.						
thickness 0.0016"						
height	SM050 3.2 mm w/ext 4.2 mm	SM100 4.6 mm	SM150 3.8 mm w/ext 5.4 mm	SM175 5.5 mm	SM200 6.4 mm	SM300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
SM050	NA	60
SM100	50	100
SM150	30	60
SM175	50	100
SM200	50	100
SM300	30	60

AU Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Composi-Tight Gold® The Gold bands are the same as the Composi-Tight® bands.						
thickness 0.0015"				NA		
height	AU050 3.2 mm w/ext 4.2 mm	AU100 4.6 mm	AU150 3.8 mm w/ext 5.4 mm	NA	AU200 6.4 mm	AU300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
AU050	25	50
AU100	50	100
AU150	25	50
AU200	50	100
AU300	25	50

M Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Composi-Tight® Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.						
thickness 0.0015"						
height	M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
M050	25	50
M100	50	100
M150	25	50
M175	50	100
M200	50	100
M300	25	50

B Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps
Composi-Tight® Original The original bands have a slightly shorter length.						
thickness 0.0015"				NA		
height	B050 3.2 mm w/ext 4.2 mm	B100 4.6 mm	B150 3.8 mm w/ext 5.4 mm	NA	B200 6.4 mm	B300 6.4 mm w/ext 8.6 mm

Item	Standard (-M)	Bulk
	Qty.	Qty.
B050	25	50
B100	50	100
B150	25	50
B200	50	100
B300	25	50

To order standard packs, add a -M after the item #.

STRATA-G™

Flexibility for a better cervical seal!

Strata-G™ Wedges utilize Soft-Face™ materials and advanced design to promote a better fit.



Stack them up.

The smooth tunnel design accommodates stacking of wedges. This feature is extremely helpful in cases of gingival recession or when simultaneous wedging from the buccal and lingual is desired.



Gentle on the papilla.

Smooth tunnel design and rounded edges glide over the papilla with less trauma.



No back-out.

Soft retentive fins smoothly fold down during wedge insertion and spring back when clear of the interproximal space, locking Strata-G™ Wedges in place.

A better fit by design.

The thoughtful design and advanced materials of Strata-G™ Wedges promote a better cervical seal.

The hollow underside sits lower and seals deeper interproximally, while the soft fins allow Strata-G™ Wedges to adapt to concavities and irregularities.



Strata-G™ Wedges

Description

Strata-G™ Wedge Kit
Strata-G™ Wedge Kit

Item

SGK4
SGK4-M

Qty.

400
200

Refills

X-Small Yellow
Small Blue
Medium Orange
Large Green

Item

SGYL-M
SGBL-M
SGOR-M
SGGR-M

Standard

Qty.

50
50
50
50

Bulk


Qty.

SGYL 100
SGBL 100
SGOR 100
SGGR 100

Helpful tip: If you have used Garrison® wedges in the past, you may find it helpful to use a size smaller than before due to the Strata-G™ flared design.

Composi-Tight® 3D Fusion™ Ultra Adaptive Wedges

Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.



Description

Composi-Tight® 3D Fusion™ Wedge Kit
Composi-Tight® 3D Fusion™ Wedge Kit

Item Qty.

FXK4 400
FXK4-M 200

Standard

Bulk

Refills

Item	Qty.
X-Small Yellow	50
Small Blue	50
Medium Orange	50
Large Green	50

Item Qty.

Item	Qty.
FXYL	100
FXBL	100
FXOR	100
FXGR	100



G-Wedge™

Perfect wedge. Perfect contact.

If you prefer traditional wedge placement with either a hemostat or cotton pliers, G-Wedges give you the same great shape as WedgeWands®, just without the applicator wand.



GSWYL

GSWBL

GSWOR

GSWGR

Description

G-Wedge™ Kit

Item Qty.

GWK4 400

Standard

Bulk

Refills



Item	Qty.
X-Small Yellow	100
Small Blue	100
Medium Orange	100
Large Green	100

Item Qty.

Item	Qty.
GWYL	300
GWBL	300
GWOR	300
GWGR	300

A+Wedge®

Controls Gingival Bleeding



GWAYL

GWABL

GWAOR

GWAGR

Description

A+Wedge® Kit

Item Qty.

GWAK4 200

Bulk


Refills

Item	Qty.
X-Small Yellow	100
Small Blue	100
Medium Orange	100
Large Green	100

Controls gingival bleeding with an aluminum sulfate coating that has an astringent effect on the gingiva, improves clarity of restoration area and eliminates need to apply an astringent separately

WedgeWands®

WedgeWands® are designed to twist off the wand, but many prefer to leave it attached to help retract the cheek and for easy removal after the procedure is complete.



WYL

WBL

WOR

WGR

Description

WedgeWands® Kit

Item Qty.

WK4 400

Standard

Bulk

Refills

Item	Qty.
X-Small Yellow	100
Small Blue	100
Medium Orange	100
Large Green	100

Item Qty.


Item	Qty.
WYL	300
WBL	300
WOR	300
WGR	300

Great technique tips and product help are only a click away!

Join thousands of subscribers on our YouTube channel and get notified instantly when new content is posted.

The Garrison® YouTube channel is NOT a collection of infomercials. It's all about helpful information on a wide variety of topics. From short two-minute tips on how to reduce flash on Class IIs, to more in depth live patient case presentations, the Garrison® channel delivers great content.

Go to youtube.com/garrisdental or search Garrison Dental on YouTube and subscribe today!



Speaking of great content...

Get CE credits direct from Garrison®!

Go to the Learning Center at garrisdental.com to see a list of pre-recorded and upcoming live webinars—many with available CE.

Register for an online account on our website to sign-up for CE courses and manage your CE credits.



Garrison Dental Solutions
Nationally Approved PACE Program Provider for FAGDIMAGD credit.
Approval does not imply acceptance by any regulatory authority or AGD endorsement.
9/1/2022 to 8/31/2026.
Provider ID# 304999

FENDERWEDGE®

Interproximal Tooth Shields

Pre-separate and Protect

FenderWedge® prevents damage to the adjacent tooth and gingiva that can occur during Class II preparations.

FenderWedge® is a combination of a steel plate and a plastic wedge. Inserted into the interproximal space, it protects the tissue and separates the teeth, thus simplifying the application of a matrix.



FenderWedge®

Description
FenderWedge® Kit

Item **Qty.**
TSK01 80

Refills	Item	Standard		Bulk	
		Qty.	Item	Qty.	Item
X-Small Purple	TSPR-M	30	TSPR	60	
Small Orange	TSOR-M	30	TSOR	60	
Medium Green	TSGR-M	30	TSGR	60	
Large Yellow	TSYL-M	30	TSYL	60	

FenderWedge is NOT a sectional matrix system nor does it replace the use of a matrix band.

Clinical Tip



Place FenderWedge® prior to preparation to pre-separate and protect the adjacent teeth.

The plastic wedge—not just the metal shield—should firmly engage the interproximal surfaces of the teeth. This will provide the optimum pre-separation force and hold the FenderWedge® securely in place (Figures 1 and 2).

Figure 3 clearly shows marks from the bur along the shield. Much better to have those on the shield than on the adjacent tooth!

Once the preparation is complete, remove the FenderWedge® and place your matrix.

Dentistry and photos courtesy of Dr. Rob Ritter, Jupiter, Florida



Visit garrisdental.com/ learning to see Dr. Ritter's full case presentation.



IPR Management in Clear Aligners Therapy

Dr. Tommaso Weinstein, DDS, PhD

Interproximal space management is a very important issue in clear aligners therapy. This case presented crowding in upper and lower arches. It was treated with clear aligner therapy through a digital planification.

The interproximal reduction (IPR) was carefully planned. The multiple thicknesses of FitStrip™ allowed IPR to be performed with predictability and control. The patient reported excellent comfort during the procedure, better than what was reported with traditional strips.



Fig 4-7. Digital planification, before (4 and 5) and after (6 and 7). Numbers refer to IPR required.



Fig 1, 2 and 3. Initial situation, upper and lower arch.



Fig 8. Use of FitStrip™ for IPR.



Fig 9. Clear aligners in place. **Fig 10 and 11.** Direct composites on incisal edges.



Fig 12 and 13. Final result in upper and lower arch.



Flexible Interproximal finishers with a Twist

FitStrips are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison® Dental Solutions**.

Do the Twist! Twist the color-coded barrel to adjust FitStrips' curvature and it's automatically locked in place.

One handy handle! FitStrip's simple, attachable handle makes a world of difference for both clinician ease-of-use and patient comfort.



Ready for a marathon!
FitStrips are very durable and long lasting. They stay sharp and ready to go even after multiple autoclave cycles.

Starts bold – finishes clean

1. Medium grit for fast removal
2. Fine grit smooths things out
3. Super fine strips put the perfect finishing touch on things

FitStrip™ Kit
Comparison
What FitStrip Kit
Should I Use?





Flexible
Interproximal finishers with a
Twist



The FitStrip™ kit that started it all!

- Primary Function: Contouring Composites and IPR up to .30mm.
- Durable diamond abrasive on just one side to protect adjacent teeth
- The FitStrip™ serrated saw is the PERFECT instrument to separate accidentally bonded teeth and reduce over-tight contacts

FitStrip™ Starter Kit

FPSK01 contains: 2 handles, 2 serrated, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm) 4 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)



Most popular!

Your Complete Solution for Interproximal Reduction (IPR)!

- Primary Function: IPR when .40mm and .50mm are needed, yet can still be used to contour composites
- IPR gauges for exact reduction
- Attachable handle for greater ease of use and patient comfort

Universal Kit

FPSK04 contains: 2 handles, 1 IPR gauge, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm) 6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)



Subgingival FitStrips are the ideal solution for the cervical area!

- Primary Function: Use with Composites, both for contouring and accessing cervical areas to remove overhangs or excess cements.
- Only 2.5mm height keeps abrasive below the contact point
- Safety zone protects contacts during insertion
- Extra thin at .05mm

FitStrip™ Subgingival Finishing and Single Sided Kit











FPSK08 contains: 2 handles, 10 Strips: 2 subgingival super fine, 2 subgingival fine, 1 subgingival medium, 1 serrated saw, 1 each single sided .08 mm, .10 mm, .13 mm, .18 mm



FitStrips – a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison Dental Solutions.**

Kits All kits come in a package of 10 FitStrips with 2 handles.

Description/Contents	Item	Primary Use
Starter Kit for finishing/polishing 2 serrated 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm) 4 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm)	FPSK01 (shown on page 18)	Finishing composites
Universal Kit for IPR and finishing/polishing 1 IPR gauge 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm) 6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)	FPSK04 (shown on page 18)	Finishing composites and IPR
IPR Kit 1 IPR gauge 10 double sided (2 of each .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)	FPSK05	IPR
Subgingival Finishing Kit 3 subgingival super fine, 3 subgingival fine, 3 subgingival medium, 1 serrated saw	FPSK07	Subgingival
Subgingival Finishing and Single Sided Kit 2 subgingival super fine, 2 subgingival fine, 1 subgingival medium, 1 serrated saw, 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)	FPSK08 (shown on page 18)	Subgingival and finishing composites

Refills		Single Sided	Double Sided
Super fine 15 Micron		FPSXFSS 0.08 mm thick	FPSXFDS 0.11 mm thick
Fine 30 Micron		FPSFSS 0.10 mm thick	FPSFDS 0.15 mm thick
Medium 46 Micron		FPSMSS 0.13 mm thick	FPSMDS 0.21 mm thick
Coarse 90 Micron		FPSCSS 0.18 mm thick	FPSCDS 0.30 mm thick
Serrated		FPSSAW 0.05 mm thick	
Medium 46 Micron			FPSTMS 0.40 mm thick
Medium 46 Micron			FPSXTMS 0.50 mm thick
Super fine 15 Micron		FPSGXFSS 0.08 mm thick	N/A
Fine 30 Micron		FPSGFSS 0.10 mm thick	N/A
Medium 46 Micron		FPSGMSS 0.13 mm thick	N/A
FitStrip™ replacement handles (2-pack) FPSHANDLE-2			
IPR Gauge FPSGAUGE			






Intra-oral Adjustment/Polishing
for Zirconia and Lithium Disilicate

Adjust and Polish Zirconia and
Lithium Disilicate the **eZr** Way!

eZr Intra-oral Adjustment:



- **No chatter, chipping or microfracture** - Unlike diamonds bonded directly to a hard metal shaft, eZr's rubber substrate cushions the diamonds to prevent damage to the restoration.
- **No heat build-up** - The carefully engineered substrate does not transfer heat during use.
- **FAST!** - eZr's coarse grit blue grinders make your adjustments fast and convenient.

eZr Intra-oral Polishing:



- **Unparalleled smooth polish** - Don't take our word for it! - "The eZr system is the only one I've tried that polishes zirconia to smooth!" Todd Brower, DDS, Holland, Michigan.
- **Increased physical stability** - Highly polished zirconia and lithium disilicate is stable and non-abrasive to opposing dentition.
- **No paste required** - Forget the mess and expense of polishing pastes. eZr™ is all you need.
- **No water required** - No need to use water when polishing. Water is only required for high speed.



1

Start by gross adjusting and contouring using the blue coarse grit grinders. Use the flat end taper for occlusal/labial/lingual/buccal surfaces and the disc for interproximal surfaces.



2

Use the green medium grit polishers for smoothing and pre-polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.



3











Use the orange fine grit polishers for high gloss polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.



Beautiful results!



Slow speed right angle latch shank
Steam autoclaveable

Item	Description	Qty.
FPZK01	Complete system kit Gross removal, smoothing, pre-polishing, and high gloss polishing	2 grinders 6 polishers
FPZC030	 Blue coarse grit disc Gross removal	5/pk
FPZC010	 Blue coarse grit flat end taper Gross removal	5/pk
FPZM010 (not in kit)	 Green medium grit small point polisher Smoothing and pre-polishing	5/pk
FPZM020	 Green medium grit large point polisher Smoothing and pre-polishing	5/pk
FPZM030	 Green medium grit disc polisher Smoothing and pre-polishing	5/pk
FPZM040	 Green medium grit cup polisher Smoothing and pre-polishing	5/pk
FPZF010 (not in kit)	 Orange fine grit small point polisher High gloss polishing	5/pk
FPZF020	 Orange fine grit large point polisher High gloss polishing	5/pk
FPZF030	 Orange fine grit disc polisher High gloss polishing	5/pk
FPZF040	 Orange fine grit cup polisher High gloss polishing	5/pk

Images courtesy of Dr. Riccardo Becciani, Italy



Direct cusp covering of vital element with a very deep cervical margin

Dr. Daniele Gensini

The case presents the need to seal the cervical margin, which is in a very deep position even if it can be isolated with a dam, and to recreate an adequate emergence profile. Thanks to the particular design of the Slick Bands™ matrices, this step was easy and made it possible to bring the restorative margin into a favorable position. Subsequently, after having recreated the direct covering of the missing cusp thanks to a previously created silicone index, the missing proximal wall was reconstructed with the aid of a Composit-Tight® 3D Fusion™ anatomical matrix, which was very useful for obtaining a correct emergence profile and an appropriate proximal contact point. Finally, the modeling of the restoration was completed and the correct occlusal contacts were checked.



Figure 1. Occlusal view of prepared tooth with cervical margin.



Figure 2. Slick Bands™ matrix placement with correct emergence profile.




Figure 3. Highly filled flowable composite for reconstruction of the deep cervical margin.



Figure 4. Restorative margin in a favorable position and with correct profile.



Figure 5. Correct adaptation of Composit-Tight® 3D Fusion™ matrix to the new restorative margin.




Figure 6. Lateral view of the correct emergence profile recreated.




Figure 7. Occlusal check.

Predictably relocate deep margins for better final restorations.

The conservative deep margin elevation procedure allows for the maximum conservation of natural tooth structure while helping to avoid more invasive procedures such as crown lengthening. Relocated supragingival margins facilitate easier scanning/impressions and proper isolation for improved bonding of inlays, onlays and direct final restorations.



Clinical Case: The margin elevation technique allows proper isolation and elevation of deep margin areas prior to restoring with an indirect partial coverage restoration (inlay or onlay). The use of the Garrison® Reel Matrix™ margin elevation band allows proper adaptation and seal of these deep areas so that a composite margin elevation can be performed.



Pre-op X-ray shows the margin location for #17 subgingivally prior to elevation.



Seat the margin elevation band subgingivally. The narrowness of the matrix will allow it to slide subgingivally and seal the margin efficiently.



Seal the dentin and elevate the deep margin with restorative composite.



Remove the matrix and finish with a #12 blade to remove any flash. The deep margin has been relocated 1-2 mm supragingival for predictable isolation and delivery of an indirect bonded restoration.



Post-op X-ray shows the completed margin elevation prior to final restoration being completed.

Photos courtesy of Matthew A. Nejad DDS, Beverly Hills, CA

ReelMatrix

Real Easy. Real Tight. Real Fast!

For all posterior restorations!

1 Real Easy:
The reels pop into the handle without any fuss! Press the reel into the handle and retract the grippers to lock the reel in place. Plus, if you use the regular matrices (yellow hubs) they can be easily pushed through contacts.

2 Real Tight:
Tight contacts and excellent contour are what Reel Matrix is all about. These thin matrices are pre-contoured in three dimensions and when combined with anatomical wedges and separator rings, the result is amazingly tight contacts at the natural height of contour.

3 Real Fast:
Pre-loaded reels mean you can just lock them in and go. The retainerless design greatly speeds up quadrant work while improving visibility and patient comfort. Add a SG500 separator ring and you virtually eliminate buccal and lingual flash.



See page 22
for Dr. Gensini's
Margin Elevation case.




RM500H
Margin elevation



RM100D
5mm deadsoft



RM200D
6.3mm deadsoft



RM100H
5mm regular



RM200H
6.3mm regular



RM100B
5mm blue



RM200B
6.3mm blue



DME-01



RMK05



Margin elevation band is .002" thin. Super-thin dead-soft bands are .0012" thin. Regular bands are .0015" thin and can be placed through existing contacts. A Blue View cure-through option at .0029" thin is available if you prefer trans-enamel curing.

Item	Description
DME-01	Deep Margin Elevation Kit 50 RM500H Margin Elevation band reels 1 RMH-100 Reel Matrix™ tensioning instrument 5 FXH150 3D Fusion™ Firm bands 10 FXH200 3D Fusion™ Firm bands 20 FXH300 3D Fusion™ Firm bands 20 FXYL 3D Fusion™ wedges 20 FXBL 3D Fusion™ wedges 20 FXOR 3D Fusion™ wedges 20 FXGR 3D Fusion™ wedges
RMK05	Kit Contains 80 matrices (20 RM500H, 10 ea. RM100D, RM200D, RM100H, RM200H, RM100B, RM200B) 1 tensioning instrument 1 PerForm™ instruments (set of 2)

Refills		Qty.
RM500H	Margin elevation band	50
RM100D	5mm dead-soft matrices	50
RM200D	6.3mm dead-soft matrices	50
RM100H	5mm regular matrices	50
RM200H	6.3mm regular matrices	50
RM100B	5mm blue matrices	50
RM200B	6.3mm blue matrices	50
RMR100	5.9mm empty reel	100
RMR200	6.7mm empty reel	100
RMH-100	Tensioning instrument	

IMPORTANT ACCESSORIES

PerForm™ Proximal Contact Instrument
(included in RMK05 kit)

Great for core build-ups

SlickBands™

Tofflemire-style

Non-stick technology developed by Garrison®: Modern bonding agents can adhere to uncoated matrix bands—Slick Bands™ slide right out.

Exclusive Garrison® shapes!



Proprietary Right Curve™ design.

The increased curvature of the **Right Curve™ matrix band** produces outstanding gingival seal and reaches over to adjacent teeth easier for broader, deeper contacts.



This creates a significantly increased funnel-shaped cone when wrapped around a tooth. It's the right choice for many restorations.

Margin Elevation Band

The deep margin elevation technique is a useful non-surgical approach to treating deep localized margins.



A margin elevation matrix will give you the best gingival adaptation and contour for the deep margin areas. It is trimmed in height to allow improved adaptation to deep margin areas.

Slick Bands™ Tofflemire-style Matrices




Clinical Tip: Use dead-soft green bands when contact is COMPLETELY broken with adjacent teeth for maximum burnishability. When contact is NOT broken with adjacent teeth, use regular gray bands.

IMPROVED Dead-soft Bands

with non-stick coating

Regular Bands

with non-stick coating



Tofflemire-style matrices are .0014" thick
Margin elevation matrices are .002" thick



See page 22
for Dr. Gensini's
Margin Elevation case.



Slick Bands™

Description

Slick Bands Kit
300 Assorted Slick Bands
100 G-Wedges
PerForm instruments (1 set; shown on page 33)

Item

SMTK06

Dead Soft Bands

Item	Standard (-M)	Bulk
	Qty.	Qty.
#1 Universal Tofflemire	50	100
Right-Curve	50	100
#13 Pedodontic	50	100
#2 Subgingival	50	100

Item

SMT200D10
SMT400D10
SMT050D10
SMT300D10

Standard (-M)	Bulk
Qty.	Qty.
50	100
50	100
50	100
50	100

Regular Bands

Item	Standard (-M)	Bulk
	Qty.	Qty.
#1 Universal Tofflemire	50	100
Right-Curve	50	100
#13 Pedodontic	50	100
#2 Subgingival	50	100

Item

SMT200H10
SMT400H10
SMT050H10
SMT300H10

Standard (-M)	Bulk
Qty.	Qty.
50	100
50	100
50	100
50	100

Regular Bands

Item	Standard (-M)	Bulk
	Qty.	Qty.
Margin elevation band	50	100

Item

SMT500H10

Standard (-M)	Bulk
Qty.	Qty.
50	100

To order standard packs, add a -M after the item #.

Cases using anteriors/cervicals



Images courtesy of Dr. Zunzaren, DDS



Images courtesy of Dr. Abdi Sameni, DDS

VariStrip™

How can you produce a natural curvature for anterior composites and eliminate black triangles? Check out Garrison's curved "Varistrips" to recreate the natural interproximal anatomy for anterior teeth.

- Helps avoid incorrect flat embrasures, black triangles and food traps that can be caused by old-fashioned straight mylars.
- No need for wedges on subgingival fillings, as the contours naturally close off the apical aspects for you.
- Excellent for class IVs, diastema closures, etc.
- Variable height allows you to slide each strip in until it precisely matches the tooth height before cutting off excess.
- The blue tint does not affect cure, but ends the annoyance of losing clear strips on patient bibs or bracket tables.

Blue View VariStrips—a must-have super-simple product for every dental office.

Dr. Michael Curtis, Bridgeport, CT, USA



Clinical Tip

After using your preferred bonding procedure, place the restorative material into the preparation in advance of using the cervical matrix.

Place the cervical edge of the matrix against the tooth (Fig. 1).

While pressing the gingival edge of the matrix firmly against the tooth, roll the matrix onto the tooth (Fig. 2). Excess material will be expressed from the occlusal edge of the matrix rather than subgingivally, greatly simplifying clean up.

Cure directly through the matrix.

For product information on the cervical kit, refer to page 32.




Photo courtesy of David S. Hornbrook, DDS, FAACD International Lecturer and Private Practice, San Diego, CA

Our Most Popular Anterior Matrix

BlueViewTM VariStripTM

Contoured Anterior Matrices



The Blue ViewTM VariStripTM is a revolutionary contoured anterior matrix providing the optimum curvature and band height for virtually every anterior restoration. The 0.0020" thin plastic anatomical strip is tapered from one end to the other. You can place the strip interproximally and then slide it until the tooth height is perfectly matched. The pre-contouring easily recreates the occluso-gingival anatomy, avoiding flat embrasures.

Curved Form/tapered: Provides the optimum curvature and band height for virtually every anterior restoration.

Blue Tint: Provides contrast between matrix and tooth structure without compromising composite polymerization.

Blue ViewTM VariStripTM

Item	Qty.
VS01	100
VS02	50



Case using the Fusion™ Anterior Matrix System

Dr. JD Corey, DDS, Muskegon, Michigan



Figure 1. Flexible matrix systems have a tendency to deform subgingivally. This can lead to non-hygienic restorations and decay on adjacent teeth (as demonstrated in this pre-operative radiograph).



Figure 2. Rubber dam is in place and teeth have been pumiced before removal of the decay.



Figure 3. Decay has been removed, the preparations have been refined, treated with gluma, etched, and bonded. These steps were completed before placing the matrix bands and wedge to ensure proper coverage of the restorative materials to the preparations.



Figure 4. Stacked matrix bands have been placed to provide facial access and are ready to be engaged by the wedge.



Figure 5. Completed restoration after removal of the rubber dam.



Figure 6. A post-operative radiograph showing a broadened interproximal contact, improved emergence profile, and appropriate marginal integrity of the restorations.

To see more cases, go to the Learning Center at garrisondental.com/learning


Fusion™ Anterior Matrix System

The **Fusion™ Anterior Matrix System** from Garrison® gives you the materials and techniques to restore many types of anterior lesions with confidence.



ANK01 Contains:

- 50 short bands
- 50 tall bands
- 25 extra small wedges
- 25 small wedges
- 25 medium wedges
- 25 large wedges



Fusion™ Anterior Bands

Refills	Item	Standard		Bulk	
		Qty.	Item	Qty.	Item
Short	AN100-M	50	AN100	100	
Tall	AN200-M	50	AN200	100	



Learn more about the Garrison® Fusion™ Anterior Matrix System



Clinical Case using the Garrison® Fusion™ Anterior Matrix System with Dr. JD Corey



Fusion™ Anterior Wedges

Description	Item	Standard		Bulk	
		Qty.	Item	Qty.	Item
Fusion Anterior Wedge Kit			ANK4	400	
Fusion Anterior Wedge Kit			ANK4-M	200	
X-Small Yellow	ANYL-M	50	ANYL	100	
Small Blue	ANBL-M	50	ANBL	100	
Medium Orange	ANOR-M	50	ANOR	100	
Large Green	ANGR-M	50	ANGR	100	

BlueView™

Cervical Matrices

"I love the Blue View Cervical Matrix! It works so well. My assistants love watching me use this system. The patients love the excitement in our faces. It makes the process easy and fast. Gotta love that! We are also hooked on using the ReelMatrix. Thanks so much for making great products."

Cynthia Brattesani, DDS,
San Francisco, CA, USA



Flexible

- adapts to individual anatomy for improved contours
- use it to directly compact the composite into the restoration
- thin edges allow excess composite to flow out for improved anatomy

Blue Tint

- provides contrast between matrix and tooth structure
- does not compromise composite polymerization

Convenient!

- protects the restoration from contamination
- eliminates time spent hand sculpting
- eliminates air-inhibited layer
- reduces finishing time
- 6 sizes for a perfect fit

KCMK2 INCLUDES

- 240 Assorted Matrices
- 1 Matrix placement Instrument

MATRICES INCLUDED:

- Anterior/Premolar
KCMA60 – 50 pcs
- Molar
KCMA70 – 50 pcs
- KCMA10 – 35 pcs
- KCMA20 – 35 pcs
- KCMA30 – 35 pcs
- KCMA40 – 35 pcs

Blue View™ Cervical Matrices

Item

KCMK2

Description

Assortment Kit w/instrument

Refills

KCMA10

Refill 100 pack

KCMA20

Refill 100 pack

KCMA30

Refill 100 pack

KCMA40

Refill 100 pack

KCMA60


Refill 100 pack

KCMA70


Refill 100 pack

KCMINST1

Placement Instrument



How PerForm™ works:




1A

STEP 1:
Place the first composite increment. Press PerForm™ into the uncured composite increment and either push (distal; fig. 1A) or pull (mesial; fig. 1B).



1B



2

STEP 2: Place the curing light directly on the Fresnel lens and cure.

PerForm

Proximal Contact Instruments

PerForm puts you in control.

- Improves contacts with any matrix system
- Control the height and breadth of contacts
- Increase the tightness of contacts
- Unique lens controls curing light diffusion
- Steam autoclave only

PerForm™

Item	Kit	Qty.
CFK01	PerForm (1/ea. size)	2/pk
CF100	1.1 mm Small Instruments	2/pk
CF200	2 mm Large Instruments	2/pk



"To ensure a tight contact a bulk fill flowable was placed in one proximal box at a time and the mylar held against the adjacent tooth with a PerForm™ proximal contact instrument.

This ensures that the contact is tight and well-formed while light curing through it."

*Dr. Jack D. Griffin, Jr., DMD, MAGD
Weldon Spring, MO*

Composi-Tight® Matrix Forceps


Precise Placement, Unbeatable Grip.



Check out the GRIP! Effortlessly remove matrix bands that otherwise would be a struggle.



Precise Placement: The tips have a 75-degree angle with an occlusal-gingival orientation, providing a more appropriate placement angle than cotton forceps. The petite, carbide tips allow the matrix to be grasped at its upper edge to prevent distortion of the interproximal area of the band.



Unbeatable Grip: The tungsten carbide tips provide exceptional grip and tactile control for tough access cases.


Matrix Forceps

BFA-1

Angled tip for occlusal placement

Garrison® Premier Ring Placement Forceps

- Extra durable
- Designed for all hand sizes and hand strengths
- Improved placement angle
- Traditional or butterfly ring placement




Garrison® Premier Ring Forceps – Improved!
FXP01




Ring Placement Forceps

Holds all your sectional matrix system rings.



Ring Placement Forceps
AUMRDF-100

Universal Forceps
for placing rubber dam clamps and rings
UNVF-100



Composite Instruments Nonstick titanium-nitride coated

Our multi-function instruments have double condenser configurations and carving blades on opposite ends. Now you can place, pack, shape and carve your way to perfectly contoured composite restorations without ever having to switch instruments!

- Large-diameter, lightweight handles for improved ergonomics.
- Provides non-stick surface for placement and shaping of composites.
- 40% smoother than stainless steel for a non-porous, non-stick surface.
- Withstands diverse cleaning/sterilization methods, including ultrasonic cleaning.



Item **Description**

TN001 **Medium Blades**
Placement of composite material
Gingival retraction to sculpt under free margins
Burnishing matrix bands



TN002 **Medium Blade with Medium Plugger**
Contouring small Class I, II, V
Packing composite material into conservative preparations
Burnishing matrix bands



TN003 **Medium Blade with Large Plugger**
Contouring large Class I, II, V
Packing composite material into large preparations
Convex plugger tip minimizes pullback of composite material



Item **Description**

TN004 **Extra Thin Blades**
Reach under margins to hone and sculpt
Sculpting interproximally
Packing cord in gingival crevice



TN005 **Long Extra Thin Blades**
Allows additional access for hard-to-reach
posterior contouring and shaping
Sculpting interproximally
Packing cord in gingival crevice











TN006 **Small/Medium
Multi-Function Plugger**
Placement of composite material in smaller Class II and III
Packing material toward gingival seats
Carving occlusal anatomy with acorn carver



Matrix Selection Guide

Being the Matrix Company, one of the most common questions we get is “which matrix do I use where?” While personal preference will drive much of the decision process, below is a chart that can provide a good starting point.

Procedure	Recommended Materials	Page	Special Notes
Distal of canine, short or malpositioned 	Strata-G™ Sectional Matrix Kit	4–5	MOST POPULAR! Separator rings designed to hold firmly on tough cases.
	Reel Matrix™—Metal or Blue View	24–25	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Slick Bands™ Tofflemire-Style Matrices	26–27	Traditional Tofflemire-style with non-stick coating.
Conservative Class II 	Strata-G™ Sectional Matrix Kit	4–5	MOST POPULAR! Superior contacts with minimal flash.
	Reel Matrix™—Metal or Blue View	24–25	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Slick Bands™ Tofflemire-Style Matrices	26–27	Traditional Tofflemire-style with non-stick coating.
Slightly wide and/or deep 	Strata-G™ Sectional Matrix Kit	4–5	MOST POPULAR! Superior contacts with minimal flash.
	Reel Matrix™—Metal or Blue View	24–25	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Slick Bands™ Tofflemire-Style Matrices	26–27	Traditional Tofflemire-style with non-stick coating.

Procedure	Recommended Materials	Page	Special Notes
Large Class II with missing cusps 	Strata-G™ Sectional Matrix Kit—Green Wide Prep Ring	4–5	EXCLUSIVE! Dentistry's first Wide Prep ring greatly simplifies these more challenging restorations.
	Reel Matrix™—Metal or Blue View	24–25	Retainerless Tofflemire-style matrices, contoured and easy to apply.
	Slick Bands™ Tofflemire-Style Matrices	26–27	Traditional Tofflemire-style with non-stick coating.
Small to medium size Class III 	Fusion™ Anterior Matrix System	30–31	Bands and wedges with ideal gingival-incisal curvature for many types of anterior lesions.
	VariStrip™ Contoured Anterior Matrices	29	Pre-contoured in three dimensions AND tapered to provide the ideal band height and contour.
Large Class III or Class IV 	Fusion™ Anterior Matrix System	30–31	Bands and wedges with ideal gingival-incisal curvature for many types of anterior lesions.
	VariStrip™ Contoured Anterior Matrices	29	Pre-contoured in three dimensions AND tapered to provide the ideal band height and contour.
Class V 	Blue View™ Cervical Matrix	32	Contours and directly compresses resin while eliminating air-inhibited layer.
Crown build-ups, post and core 	Reel Matrix™	24–25	Best choice. Margin Elevation band available for very deep buildups. Retainerless convenience.
	Slick Bands™ Tofflemire-Style Non-Stick Matrices	26–27	Margin Elevation band available for very deep buildups if you prefer a traditional matrix retainer.

Customer Service Information

Call **+61 416 232 252**

Email gds@garrisondental.com

On-line **garrisondental.com**
Look for our Learning Center
with clinical tips, case studies
and technique guides.

Mail 150 DeWitt Lane
Spring Lake, MI 49456
USA

Distribution

ivoclar

NEW ZEALAND
www.ivoclar.com

Orders: 0508 IVOCLAR
(486 252)

To find a dealer outside of Australia or New Zealand contact our headquarters in
Spring Lake, MI USA. | Phone: +1 616.842.2244 | Email: gds@garrisondental.com



Some products may be special order only or not available through all dealers.

Garrison® Guarantee

You can buy with complete confidence when you select Garrison products. We support all our products with expert customer service, outstanding warranties, and unrivaled satisfaction guarantees. Please visit our website for details.

Garrison

150 DeWitt Lane | Spring Lake, MI 49456 USA
toll-free 888.437.0032 | www.garrisondental.com

INAUCAT24C