1 Intended use

Intended purpose

Etching of the tooth structure.

Patient target group

- Patients with permanent teeth
- Patients with deciduous teeth

Intended users / Special training

- Dentists (clinical procedure)
- Dental assistants (preparation and re-filling of the syringe)
- No special training required

Use

For dental use only.

Description

For use in the enamel-etch technique or total-etch technique (enamel and dentin conditioning).

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Areas of application

- For use in the enamel-etch technique or total-etch technique (enamel and dentin conditioning).

Contraindications

- Do not apply to areas close to the pulp.
- The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

Limitations of use

If the stipulated working technique cannot be applied.

Side effects

There are no known side effects to date.

Interactions

The surfaces of glass ionomer-based bases/liners are dissolved by Total Etch.

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics

Composition

Water, phosphoric acid, polyethylene glycol, silicon dioxide

Total Etch

[sl]	Navodila za uporabo Gel za jekkanje; 37% fosforna kislina	[sq]	Udhëzime përdorimi Xhel ashpërues; 37% acid fosforik
[hr]	Upute za uporabu Gel za jekkanje; 37% fosforna kiselina	[ro]	Instrucțiuni de utilizare Gel pentru gravaj acid; acid fosforic 37%
[cs]	Návod k použití Leptací gel; 37 % kyselina fosforečná	[uk]	Інструкція щодо використання Травильний гель; 37% фосфорної кислоти
[sk]	Návod na používanie Leptavý gel; 37% kyselina fosforečná	[et]	Kasutamisi juhend Sõõritav geel; 37% fosforhape
[hu]	Használati utasítás Marógél; 37% fosforosav	[lv]	Lietošanas instrukcija Oforta gēls; 37% fosforskābe
[sr]	Упутство за употребу Жутство за употребу Гел за нагризање; 37% фосфорна киселина	[it]	Naudojimo instrukcija Ēsdinimo gelis; 37 % fosforo rūgštis
[mk]	Упатство за употреба Гел за нагризување; 37% фосфорна киселина		
[bg]	Инструкции за употреба Етчащ гел; 37% фосфорна киселина		

Total Etch

[en]	Instructions for Use Etching Gel; 37% phosphoric acid	[fi]	Käyttöohjeet Etsausgeeli; 37% fosforihappo
[de]	Gebrauchsinformation Ätzelgel; 37% Phosphorsäure	[no]	Bruksanvisning Etsegel; 37% fosforsyre
[fr]	Mode d'emploi Gel de mordantage; 37% d'acide phosphorique	[nl]	Gebruiksaanwijzing Etsgel; 37% fosforzuur
[it]	Istruzioni d'uso Gel mordenzante; 37% acido fosforico	[el]	Οδηγίες Χρήσεως Απορροφητικό ζελέ; 37% φωσφορικό οξύ
[es]	Instrucciones de uso Gel de grabado; Acido fosfórico al 37%	[tr]	Kullanma Talimatı Asit jel; %37 Fosforik asit
[pt]	Instruções de Uso Gel para condicionamento ácido; ácido fosfórico a 37 %	[ru]	Инструкция по применению Протравливающий гель; 37% фосфорной кислоты
[sv]	Bruksanvisning Etsgel; 37%-ig fosforsyra	[pl]	Instrukcja stosowania Żel do wytrawiania tkanek zębnych; 37% kwas fosforowy
[da]	Brugsanvisning Ætsgel; 37% fosforsyre		

2 Application

- Use a suitable hygienic protective sleeve for the syringe as indicated by the manufacturer.
- Apply Total Etch on enamel and subsequently on dentin and allow a reaction time of 15 seconds.
- Afterwards, thoroughly rinse off the phosphoric acid with water spray and dry the tooth surfaces with oil-free air. Avoid excessive drying of the dentin.
- The etched enamel surface should have a chalky white appearance. If this is not the case, repeat etching with a shorter reaction time.
- If the tooth surface becomes accidentally contaminated (e.g. with saliva), repeat the etching procedure with a reaction time of max. 10 seconds.
- The reaction time on unprepared enamel (e.g. fissure sealing) is 30-60 seconds. If the enamel is selectively etched, Total Etch has to be left to react for 15-30 seconds.
- Do not apply the product to areas close to the pulp.
- **Adequate relative or absolute isolation is required.**

Application notes

- Avoid any contact of the reusable syringe not covered with a protective sleeve with the mouth of the patient.
- Dispose of contaminated syringes.

Refilling empty Total Etch application syringes

It is necessary to wear gloves and protective goggles for safety reasons.

- Empty 2-ml syringes can be conveniently refilled with the contents of the Total Etch jumbo syringe up to 15 times.
 - Remove the screw-on cap from the 30-gram jumbo syringe and the empty 2-ml syringe.
 - Join the tips of the two syringes using a twisting, clockwise motion. Make sure that the two syringes are firmly connected.
- Caution:** Do not strip the thread.
- Refill the 2-ml syringe by slowly and steadily pulling back the plunger of the 2-ml syringe. Simultaneously exert light pressure on the plunger of the jumbo syringe using the other hand to support the refilling process.
 - Disconnect the two syringes using a counter-clockwise twist.
 - Replace the screw-on cap on the 30-gram jumbo syringe.
 - In order to ensure traceability of batches, write the LOT number and date of expiration of the 30-gram jumbo syringe on one of the labels included in the delivery form and stick it on the filled 2-ml syringe.
 - Attach an application tip to the 2-ml syringe before using it.

3 Safety information

- Avoid overfilling the 2-ml syringe.
- Under no circumstances should you exceed the plunger safety stop of the small syringe due to the risk of acid burn!
- Do not point the syringes at people during refilling.
- Do not fill etchant from the 2-ml syringe back into the jumbo syringe.
- Use only the Total Etch application tips delivered with this package. The use of other application tips may have an adverse effect on the viscosity of the etchant.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: www.ivoclar.com, and your responsible competent authority.
- The current Instructions for Use are available in the download section of the Ivoclar Vivadent AG website (www.ivoclar.com).

Warnings

Observe the Safety Data Sheet (SDS) (available in the download section of the Ivoclar Vivadent AG website www.ivoclar.com).

- Phosphoric acid is corrosive.
- Avoid contact with the eyes, mucous membrane and skin (goggles are recommended for both patients and operators).
- If the material should accidentally come in contact with the eyes, immediately rinse with copious amounts of water and consult a physician / ophthalmologist.
- After skin contact, immediately rinse with water.

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks. Some of these risks are listed below:

- Postoperative sensitivity due to over-etching of the dentin

4 Shelf life and storage

- Storage temperature 2–28 °C (36–82 °F)
- If the indicated storage temperature is exceeded, proper functioning of the product cannot be ensured.
- Do not use the product after the indicated date of expiration.
 - Expiry date: See note on the packages.

5 Additional information

Keep material out of the reach of children!

Not all products are available in all countries.

The product has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of use. The user is responsible for testing the product for its suitability and use for any purpose not explicitly stated in the Instructions.

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www.ivoclar.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

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www.ivoclar.us

Ivoclar Vivadent Shanghai Trading Co., Ltd.

2/F Building 1, 881 Wuding Road, Jing An District, 200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

www.ivoclar.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 3399, Fax +57 1 633 1663

www.ivoclar.co

Ivoclar Vivadent SAS

B.P. 118, 74410 Saint-Jorioz, France

Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

www.ivoclar.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, 73479 Ellwangen, Jagst, Germany

Tel. +49 7961 889 0, Fax +49 7961 6326

www.ivoclar.de

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate

Veera Desai Road, Andheri (West), Mumbai, 400 053, India

Tel. +91 22 2673 0302, Fax +91 22 2673 0301

www.ivoclar.in

Ivoclar Vivadent Marketing Ltd.

The Icon, Horizon Broadway BSD, Block M5 No. 1

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www.ivoclar.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale, PO Box 303011 North Harbour

Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

www.ivoclar.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawła II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

www.ivoclar.pl

Ivoclar Vivadent LLC

Prospekt Andropova 18 korp. 6/, office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

www.ivoclar.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204

P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 11 293 8345, Fax +966 11 293 8344

www.ivoclar.com

Ivoclar Vivadent S.L.U.

Carretera de Fuencarral nº24, Portal 1 – Planta Baja
28108-Alcobendas (Madrid), Spain
Tel. +34 91 375 78 20, Fax +34 91 375 78 38
www.ivoclar.es

Ivoclar Vivadent AB

Gustav III:s Boulevard 50, 169 74 Solna, Sweden
Tel. +46 8 514 939 30, Fax +46 8 514 939 40
www.ivoclar.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi, Sakayik Sokak, Nisantasi' Plaza No:38/2
Kat:5 Daire:24, 34021 Sisli – Istanbul, Turkey
Tel. +90 212 343 0802, Fax +90 212 343 0842
www.ivoclar.com

Ivoclar Vivadent Limited

Compass Building, Feldspar Close, Warrens Business Park
Enderby, Leicester LE19 4SD, United Kingdom
Tel. +44 116 284 7880, Fax +44 116 284 7881
www.ivoclar.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA
Tel. +1 800 533 6825, Fax +1 716 691 2285
www.ivoclar.us