Telio® CAD

[en] Instructions for Use

PMMA disc / PMMA block

[de] Gebrauchsinformation

PMMA-Scheibe / PMMA-Block

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Disque PMMA / Bloc PMMA

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Disco in PMMA / Blocchetto in PMMA

[es] Instrucciones de uso

Disco de PMMA / Bloque de PMMA

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Disco de PMMA / Bloco de PMMA

[sv] Bruksanvisning

PMMA-disk / PMMA block

[da] Brugsanvisning

PMMA-disk / PMMA blok

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PMMA-kiekko / PMMA-lohko

[no] Bruksanvisning

PMMA-disk / PMMA-blokk

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PMMA-schijf / PMMA blok

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1 Intended use

Intended purpose

Fabrication of temporary crowns and bridges

Patient target group

- Patients with permanent teeth

Intended users

- Dentists (chairside fabrication of restorations, clinical procedure)
- Dental laboratory technicians (fabrication of restorations in the dental laboratory)

Special training

None

Use

For dental use only.

Description

Cross-linked PMMA blocks and discs for the fabrication of long-term temporary crowns and bridges with up to two connected pontics in the anterior and posterior region.

Technical specifications

Property	Specifications	Typical mean value (according to ISO 10477:2020)
Flexural strength	≥ 100 MPa	130 MPa
Water absorption	≤ 40 $\mu\text{g}/\text{mm}^3$	28 $\mu\text{g}/\text{mm}^3$
Water solubility	≤ 7.5 $\mu\text{g}/\text{mm}^3$	0.6 $\mu\text{g}/\text{mm}^3$

Indications

- Missing tooth structure in anterior and posterior teeth
- Partial edentulism in the anterior and posterior region

Types of restorations:

- Temporary crowns and bridges

Contraindications

If the patient is known to be allergic to any of the product's ingredients.

Limitations of use

- Bridge constructions with more than two connected bridge pontics
- Hybrid abutments and hybrid abutment crowns
- Untreated bruxism
- Reuse of the final restoration
- Processing the product in a non-compatible CAD/CAM system
- Failure to comply with the required minimum thicknesses
- Permanent restorations

Side effects

None known to date

Interactions

None known to date

Clinical benefit

- Reconstruction of chewing function
- Restoration of esthetics
- Protection of dentin lesions

Composition

Polymethyl methacrylate (PMMA), pigments

2 Usage

NOTICE! All working steps in the oral cavity of the patient must exclusively be performed by a qualified professional (e.g. dentist).

2.1 Fabricating the restoration

NOTICE! Observe the following parameters:

Table 1

Types of restorations	Anterior region		Posterior region	
	Minimum wall thickness of the restoration	Connector dimensions	Minimum wall thickness of the restoration	Connector dimensions
Crown	occlusal 1.5 mm circular 0.8 mm	–	occlusal 1.5 mm circular 0.8 mm	–
3-unit bridge	occlusal 1.5 mm circular 0.8 mm	min. 12 mm ²	occlusal 1.5 mm circular 0.8 mm	min. 12 mm ²
4- and multi-unit bridges with 2 connected pontics	occlusal 1.5 mm circular 0.8 mm	min. 12 mm ²	occlusal 1.5 mm circular 0.8 mm	min. 16 mm ²

Determining the tooth shade

1. Clean the teeth professionally.
2. Determine the tooth shade with the tooth still moist using the Vita A-D shade guide.

Preparing the tooth according to established preparation guidelines

1. Prepare tooth while adhering to the minimum wall thicknesses, see Table 1.
2. Round out inner edges.

Digitizing the oral situation

1. Scan model/oral situation intraorally.

Designing and nesting the restoration in the CAD/CAM software

1. Design the restoration while adhering to the parameters, see Table 1.
2. Nest the restoration.

Fabricating the restoration

1. Fix the product in the holder.
2. Fabricate the restoration.

Finishing the restoration

1. Remove the product from the holder.
2. Separate the attachment points.
3. Smooth out the attachment point.
4. Dentist: Optional: Perform intraoral try-in. See "Performing intraoral try-in".
5. Dental technician: Check fit of the restoration on the die.
6. Process and polish the restoration with light pressure and low speed while adhering to the parameters, see Table 1.
7. Clean the restoration.
8. Dry the restoration with oil-free compressed air.

Performing intraoral try-in

1. Check static and dynamic occlusion.
2. Clean the restoration.
3. Dry the restoration with oil-free compressed air.

OPTIONAL: Characterizing the restoration

1. Perform cut-back.
2. Sandblast the area to be built up with Al₂O₃, 80 – 100 µm, 1 – 2 bar (15 – 29 psi) or roughen it with a coarse diamond bur.
3. Clean the restoration.
4. Dry the restoration with oil-free compressed air.
5. Apply a light-curing bonding agent to the area to be built up. Observe the Instructions for Use.
6. Apply a light-curing composite material to the restoration and polymerize. Observe the Instructions for Use.
7. Process and polish the restoration with light pressure and low speed while adhering to the parameters, see Table 1.

Preparing the restoration for cementation

1. Blast the inner aspects of the restoration with Al₂O₃, 80 – 100 µm, 1 – 2 bar (15 – 29 psi).
2. Clean the restoration.
3. Dry the restoration with oil-free compressed air.
4. For adhesive cementation: Apply a light-curing bonding agent. Observe the Instructions for Use.

Cementing the restoration

NOTICE! Check long-term temporaries with temporary cementation regularly and re-cement if necessary.

NOTICE! If adhesively cemented, the restoration may be difficult to remove.

1. Temporary cementation: Cement the restoration with temporary, eugenol-free luting material.
2. Adhesive long-term cementation: Cement the restoration adhesively.

OPTIONAL: Polishing the restoration

1. Check static and dynamic occlusion.
2. Polish the restoration.

2.2 Relining the restoration

1. Sandblast the area to be built up with Al_2O_3 , 80 – 100 μm , 1 – 2 bar (15 – 29 psi) or roughen it with a coarse diamond bur.
2. Clean the restoration.
3. Dry the restoration with oil-free compressed air.
4. Apply a light-curing bonding agent to the area to be built up. Observe the Instructions for Use.
5. Fill the inner aspect with a flowable composite material. Observe the Instructions for Use.
6. Isolate the abutment with vaseline.
7. Place the restoration onto the abutment.
8. Wait for the curing time to pass. Observe the instructions of the manufacturer.
9. Remove the restoration from the oral cavity.
10. Allow the restoration to cure.
11. Remove excess material.
12. Polish the restoration.
13. Prepare the restoration for cementation, see "Preparing the restoration for cementation".
14. Cement the restoration, see "Cementing the restoration".

2.3 Adjusting/repairing the restoration

1. Sandblast the area to be built up with Al_2O_3 , 80 – 100 μm , 1 – 2 bar (15 – 29 psi) or roughen it with a coarse diamond bur.
2. Clean the restoration.
3. Dry the restoration with oil-free compressed air.
4. Apply a light-curing bonding agent to the area to be built up. Observe the Instructions for Use.
5. Apply a light-curing composite material to the restoration and polymerize. Observe the Instructions for Use.
6. Polish the restoration.
7. Prepare the restoration for cementation, see "Preparing the restoration for cementation".
8. Cement the restoration, see "Cementing the restoration".

3 Safety information

- The product has been developed solely for use in dentistry. Process according to the Instructions for Use.
- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, www.ivoclar.com and your responsible competent authority.
- During finishing, grinding dust may be generated. Grinding dust can be harmful to health. Do not inhale grinding dust. Use an extraction unit and wear protective goggles and a face mask. Wear protective gloves.
- Check to make sure that the packaging and the product are intact and undamaged before use. If in doubt, contact Ivoclar Vivadent AG or your local sales partner.

Supporting documents

Document	To be found at:
Current version of Instructions for Use	www.ivoclar.com/eIFU
Structure of Instructions for Use and Warnings	www.ivoclar.com/eIFU
Explanation of symbols	www.ivoclar.com/eIFU
Safety Data Sheet (SDS)	www.ivoclar.com

Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following known residual clinical risks exist:

- Chipping / fracture / debonding of the restorative material may lead to accidental swallowing or inhalation and dental retreatment
- Cement excess may lead to irritation of the soft tissue/gingiva. Progressive inflammation may lead to bone resorption or peri-implant disease.

3.1 Shelf life and storage

- Store the product free of vibration.
- Store the product in a dry place.
- Protect the product from sunlight.

4 Additional information

Keep material out of the reach of children!

The manufacturer accepts no liability for damages resulting from improper use of the product. Furthermore, the user is responsible for determining the product's suitability for the user's own purposes before use.

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Polymethylmethacrylat (PMMA), Pigmente

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