

PrograScan[®] PS7



Simultaneous high speed for perfection

The high-end lab scanner PrograScan PS7 scans upper and lower jaw models in one single process: simultaneously and fully automatically.

Technical Data

Camera	8 (8 MP each)
Scan tables	2
Auto-Start	Yes
Accuracy (ISO 12836)	5 µm
Repeatability	<2 µm
Texture	Colour
Scan strategy	All-in-one /die in model
Output format	STL, PLY
Voltage range	100 – 240 VAC
Frequency range	50 – 60 Hz
Rated power (max.)	90 Wh
Scanner voltage	24V DC ±10%
Type of USB port	USB 3.0 TYP C (5 Gbit/s)
Number of USB ports	2
Weight	13,6 kg
Dimensions (W x H x D)	396.0 x 436.5 x 380.5 mm

Date information prepared: 2023-11-06