

Simply good

Ivostar®

Denture teeth made of quality PMMA



Tried and tested quality and reliability

At Ivoclar we are committed to providing you with products of the highest quality standards.

Before we supply our denture teeth to you, we check them again and again until we are 200 percent satisfied with their quality. We thoroughly test our products to make sure that you can use them successfully to create highly esthetic and functional removable dentures.



The fascination of prosthetics – more than a denture

We pay attention to details, but never lose sight of the bigger picture. Our prefabricated denture teeth are an integral component of our functional and coordinated prosthodontic workflow.

Individualized and efficient – for you and your clients.



Proven material

- ✓ Cross-linked PMMA
- ✓ Outstanding chemical bond to denture base materials
- ✓ Versatile uses



200%
tested
quality



Large shade selection

- ✓ A-D, Bleach
- ✓ High level of flexibility in routine laboratory work
- ✓ Excellent performance in complete, partial and implant prosthetics



Lifelike

- ✓ Three-layer anterior tooth
- ✓ Shape-supporting individualized layering
- ✓ Clinically proven material




Your benefits:

- ✓ Rigorous quality assurance during production
- ✓ Stringent quality checking of the finished product
- ✓ 2 x 100 percent quality control for your peace of mind

At a glance	Ivostar®	Gnathostar®
Tooth moulds	21 different anterior tooth sets (15 upper jaw, 6 lower jaw)	Posterior teeth for conventional tooth-to-two-teeth set-up (5 upper jaw, 5 lower jaw)
Tooth colors	16 A-D shades 4 Bleach shades	
Material	Cross-linked PMMA	
Number of layers	3	2

Conventional removable dentures restoring individuality and function



Preliminary bite registration

Precision from the start with Centric Tray



Functional impression

Centricity to the point with Gnathometer M



Registration

Successful results in no time with the UTS 3D



Tooth set-up

Customized dentures with our prefabricated teeth



Try-in

Predictable results, verified on the patient



Finalization

Versatile denture base materials for the pouring, condensing and injection techniques



Protect

Effective preventive care with the Cervitec line of products

United States

Ivoclar | 175 Pineview Drive | Amherst | New York 14228
716-691-0010 | 800-533-6825 | Fax: 716-691-2285

Canada

Ivoclar | 1 - 6600 Dixie Road | Mississauga | Ontario L5T 2Y2
905-670-8499 | 800-263-8182 | Fax: 905-670-3102

ivoclar.com