


CAD/CAM

Block Overview

for PrograMill®




IPS e.max[®] CAD

Lithium disilicate glass-ceramic (LS₂) | Lithium-Disilikat Glaskeramik (LS₂)


HT (High Translucency)

Size | Grösse



I12 C14

Shades | Farben



A-D Bleach

A1 A1 A3 A3.5 A4 B1 B2 B3 B4 C1 C2 C3 C4 D2 D3 D4 BL 1 BL 2 BL 3 BL 4

MT (Medium Translucency)

Size | Grösse



C14

Shades | Farben




A-D Bleach

A1 A2 A3 B1 BL 2 BL 3 BL 4


LT (Low Translucency)

Size | Grösse



C14


Shades | Farben



A-D Bleach


A1 A2 A3 A3.5 A4 B1 B2 B3 B4 C1 C2 C3 C4 D2 D3 D4 BL 1 BL 2 BL 3 BL 4

Size | Grösse



C16 B32

Shades | Farben



A-D Bleach

A1 A2 A3 A3.5 B1 B2 C1 C2 D2 BL 2

MO (Medium Opacity)

Size | Grösse



C14


Shades | Farben



MO 0 MO 1 MO 2 MO 3 MO 4


Impulse

Size | Grösse



C14

Shades | Farben



Opal 1 Opal 2

Types of restorations | Restaurationsarten

Thin veneers (0,4 mm), veneers, occlusal veneers, inlays, onlays, partial crowns, crowns, minimal invasive crowns (≥1 mm), 3-unit-bridges (up to the second premolar as the terminal abutment)

Dünne Veneers (0,4 mm), Veneers, Okklusale Veneers, Inlays, Onlays, Teilkronen, Einzelzahnkronen, minimalinvasive Kronen (≥1 mm), 3-gliedrige Brücken (bis zum zweiten Prämolaren als endständigen Pfeiler)

IPS Empress[®] CAD

Leucite glass-ceramic | Leuzit-Glaskeramik

Multi

Size | Grösse



C14 L

Shades | Farben



A-D Bleach

A1 A2 A3 A3.5 B1 BL 1 BL 3


HT (High Translucency)

Size | Grösse



I12

Shades | Farben



A-D

A1 A2 A3 A3.5 B1 B2 B3 C2 D3

LT (Low Translucency)

Size | Grösse



C14

Shades | Farben



A-D Bleach

A1 A2 A3 A3.5 B1 B2 B3 C2 D3 BL 1 BL 2 BL 3 BL 4

Types of restorations | Restaurationsarten

Veneers, inlays, onlays, partial crowns, crowns

Veneers, Inlays, Onlays, Teilkronen, Kronen

Tetric[®] CAD

Composite

HT (High Translucency)

Size | Grösse



C 14

Shades | Farben



A-D Bleach

A1 A2 A3 A3.5

MT (Medium Translucency)

Size | Grösse



C 14

Shades | Farben



A-D Bleach

A1 A2 A3 A3.5 BL


Types of restorations | Restaurationsarten

Veneers, inlays, onlays (e.g. partial crowns, occlusal veneers), crowns

Veneers, Inlays, Onlays (z.B. Teilkronen, okklusale Veneers), Kronen

Block Dimensions | Block-Dimensionen

Abbreviation Abkürzung	Block	Size Grösse in mm		
		Width Breite	Length Länge	Height Höhe
IPS e.max® CAD				
Inlay	I 12	10.4	12.5	15.0
Crown/Abutment	C 14	12.4	14.5	18.0
Crown/Abutment	C 16	17.8	15.8	18.0
Bridge	B 32	14.5	14.5	32.0
IPS Empress® CAD				
Inlay	I 12	10.0	12.0	15.0
Crown	C 14	12.0	14.0	18.0
Crown	C 14 L	14.4	14.4	18.0
Tetric® CAD				
Crown	C 14	11.8	13.8	17.8



Date information prepared: 2022-07, Rev. 3
695567/EN-DE

Ivoclar Vivadent AG | Bedererstrasse 2 | 9494 Schaan | Liechtenstein
Tel. +423 235 35 35 | Fax +423 235 33 60

ivoclar.com