

Time Tested. Time Proven.
Time you tried it!



Heliomolar[®] Microfill Restoratives

Heliomolar®

The Clinically Tried-and-Tested Universal Composite for Anterior and Posterior Restorations



20 Years of Proven Performance

The reinforced microfiller formula of Heliomolar provides excellent esthetics and polishability, high strengths, and renowned wear resistance. These features, combined with a proven performance record, have made Heliomolar the benchmark composite material by which all new materials are measured.


Indicated for both anterior and posterior applications, Heliomolar provides the following features and benefits:

- Outstanding Microfill Esthetics
- Impressively Low Wear Abrasion
- Continuous Fluoride Release
- Universal Application (Posterior & Anterior)
- Syringe and Cavifil Delivery

Heliomolar is available in 11 Esthetic Shades including a pedo shade and two enamel shades.



An independent study* demonstrated the outstanding wear performance of Heliomolar.



This 15 year old Heliomolar restoration demonstrates the TIME TESTED success regarding function *and* esthetics. (Dentistry by Doug Lambert, DDS)

Heliomolar® Flow

The Flowable Microfill Developed on 20 Years of Proven Performance and Success.



The Flowable Microfill That Gets Down into the Tightest Spots.

Based on a formulation with 20 years of proven performance, Heliomolar Flow provides excellent esthetics, superior polishability and excellent flowability. Packaged in Luer-Lock syringes, each Heliomolar Flow refill includes small applicator tips (0.9mm), pre-bent at 45° to facilitate easy and precise placement.

**Indicated for all flowable applications,
Heliomolar Flow features:**

- Outstanding Microfill Esthetics
- Impressively Low Wear Abrasion
- Excellent Flow Properties
- Continuous Fluoride Release
- Optimal Wetting
- Luer-Lock Syringe Delivery

Heliomolar Flow is available in 7 Esthetic Shades including a pedo shade and two enamel shades.



Heliomolar® Flow



Heliomolar® Flow's unique wetting properties allow for complete penetration to the depth of the cavity prep.

Other Flowables



Other hybrid flowables do not routinely penetrate to the base of the cavity preparation, leaving undesirable gaps and spaces.

Heliomolar® HB

The First “Heavy Body” Reinforced Microfill



Heliomolar HB is CLEARLY the PERFECT material for posterior restorations.

Along with excellent esthetics and superior polishability, Heliomolar HB’s advanced “Heavy Body” handling characteristics provides the packability and sculptability necessary for posterior restorations.

Heliomolar HB’s two NEW translucent shades, Translucent White (TW) and Translucent Gray (TG) take esthetics to a whole new level. These virtually clear shades, combined with the exceptionally low wear rates of Heliomolar HB, make them the perfect choice for the final occlusal layer on all posterior restorations. Now, dental professionals can create posterior restorations with highly translucent esthetics and long lasting performance.




Indicated for all posterior restorations,

Heliomolar HB Features:

- “Heavy Body” Consistency
- Outstanding Microfill Esthetics
- Impressively Low Wear Abrasion
- High levels of Translucency
- 20 Second Cure Time

Heliomolar HB is available in 11 esthetic shades.



BEFORE

AFTER

Heliomolar HB restoration using Heliomolar Flow as the initial liner, Heliomolar HB in the body of the restoration and the new translucent shades of Heliomolar HB as the final occlusal layer.

(Denistry by Robert Ritter, DDS)

✓ Time Tested.

Heliomolar[®], with over 20 years of clinical success, is the material of choice for class II restorations according to a leading independent research organization.*

✓ Time Proven.

Heliomolar[®] is the benchmark composite, offering exceptionally low wear rates similar to enamel.*

✓ Time you tried it.






The Heliomolar[®] family of products provides dental professionals with a proven formulation in three viscosities satisfying all indications and placement techniques. Numerous independent studies have been published supporting the clinical success of Heliomolar's reinforced microfill chemistry over the past 20 years.

The Heliomolar[®] Family of Products

1 Proven Microfill Chemistry – 3 Viscosities

Heliomolar[®]Flow, Heliomolar[®] & Heliomolar[®]HB all feature:

- Outstanding Microfill Esthetics in a Variety of Shades
- Superior Polishability
- Excellent Handling Properties
- Renowned Wear Resistance
- High Radiopacity*
- Fluoride Release
- Multiple Delivery Forms

	Heliomolar[®] Flow Flowable Reinforced Microfill Restorative	Heliomolar[®] Universal Reinforced Microfill Restorative	Heliomolar[®] HB Heavy-Body Reinforced Microfill Restorative
Viscosity	Light Body	Regular Body	Heavy Body
Indication	<ul style="list-style-type: none"> • Small size cavities • Cervical restorations • Anterior restorations • Preventive resin restorations • Lining 	<ul style="list-style-type: none"> • Class I • Class II • Class III • Class IV • Class V 	<ul style="list-style-type: none"> • Class I • Class II
Independent Evaluation Results	 		 
*Heliomolar HB TW & TG shades are not radiopaque			

The Heliomolar® Family of Products

1 Proven Microfill Chemistry – 3 Viscosities

Heliomolar® Flow

Product Code	Product	Contents
557019	Heliomolar Flow Assortment	7 Syringes (1.6 gm each) Heliomolar Flow (1 each in shade 105, 110, 140, 210, 410, 110T, 420T) 40 Luer-Lock Applicator Tips (0.9mm)
	Heliomolar Flow	1 syringe (1.6gm)

Heliomolar® Flow

Product Code	Product	Content
	Syringe Refills	5 Luer-Lock Applicator Tips (0.9mm)
557029	Shade 105--/18	
557030	Shade 110/A1/20	
557031	Shade 140/A2/22	
557032	Shade 210/A3/U	
557033	Shade 410/D3/35	
557034	Shade 110T/A1T/20T	
557035	Shade 420T/--/38T	
556744	Heliomolar Flow Luer Lock Applicator Tips	20 Heliomolar Flow Luer-Lock Applicator Tips (0.9 mm)

Heliomolar®

Product Code	Product	Content
550576	Heliomolar 8 Syringes Assortment	(3gm each) Heliomolar (1 each shade 105, 110, 140, 210, 410, 510, 110T, 420T) 1 Heliomolar Shade Guide
	Heliomolar Syringe Refills	1 syringe (3gm) Heliomolar
550555	Shade 105/18	
532686	Shade 110/A1/20	
532687	Shade 140/A2/22	
532685	Shade 210/A3/Univ.	
532688	Shade 340/A4/24	
532689	Shade 310/B3/30	
532690	Shade 440/D4/32	
532691	Shade 410/D3/35	
532692	Shade 510/C3/36	
550556	Shade 110T/WE/20-T	
550557	Shade 420T/GE/38-T	
550575	Heliomolar Cavifil Assortment	40 Cavifils (0.25 gm each) Heliomolar (5 each shade 105, 110, 140, 210, 410, 510, 110T, 420T) 1 Heliomolar Shade Guide
	Heliomolar Cavifil Refills	20 Cavifils (0.25 gm each) Heliomolar
550558	Shade 105/18	
541501	Shade 110/A1/20	
541502	Shade 140/A2/22	
541503	Shade 210/A3/Univ.	
541505	Shade 340/A4/24	
541504	Shade 310/B3/30	
541507	Shade 440/D4/32	
541506	Shade 410/D3/35	
541508	Shade 510/C3/36	
550559	Shade 110T/WE/20-T	
550560	Shade 420T/GE/38-T	

Note: Heliomolar HB TW and TG shades are not radiopaque.

Heliomolar® HB

Product Code	Product	Content
	Heliomolar HB Refills	1 Syringe (3 gm) Heliomolar HB
559644	Shade 110/A1	
560195	Shade 140/A2	
560196	Shade 210/A3	
560197	Shade ---/B1	
560198	Shade 130/B2	
560199	Shade 430/C2	
560200	Shade 1---/D2	
560201	Shade 110T/A1T/WE	
560202	Shade 420T/--/GE	
563478	Shade TW/Translucent White	
563479	Shade TG/Translucent Gray	
	Heliomolar HB Cavifil Refills	20 Cavifils (0.25 gm each) Heliomolar HB
560066	Shade 110/A1	
560067	Shade 140/A2	
560068	Shade 210/A3	
560069	Shade ---/B1	
560070	Shade 130/B2	
560071	Shade 430/C2	
560072	Shade 1---/D2	
560073	Shade 110T/A1T/WE	
560074	Shade 420T/--/GE	
563480	Shade TW/Translucent White	
563481	Shade TG/Translucent Gray	
560203	Heliomolar HB Shade Guide	



100% CUSTOMER SATISFACTION
GUARANTEED!

ivoclarvivadent.com

© 2010 Ivoclar Vivadent, Inc. Heliomolar is a registered trademark of Ivoclar Vivadent, Inc. *Independent non-profit dental products testing institution. Data available upon request. Herculite XRV is not a registered trademark of Ivoclar Vivadent, Inc. Rev. 3 03/2010 6910092

UNITED STATES
175 Pineview Drive
Amherst, New York 14228
716-691-0010
Fax: 716-691-2285
www.ivoclarvivadent.com

CANADA
2785 Skymark Ave., Unit 1
Mississauga, Ontario L4W 4Y3
800-263-8182
Fax: 905-238-5711
www.ivoclarvivadent.com

Heliomolar®
UNIVERSAL MICROFILL COMPOSITE

ivoclar vivadent
passion vision innovation