

SAFETY DATA SHEET (GHS)

Issue Date: 30.11.2018
Version: 1.0.0
Revision Date: 28.11.2023
Print Date: 30.11.2018

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE/PRODUCT AND MANUFACTURER/IMPORTER

1.1 Product identifier:-

Product name: | IPS InLine PoM Touch Up
Product number: | Void

1.2 Other means of identification:-

Not applicable.

1.3 Recommended use of the chemical and restrictions on use:-

Not applicable.

Identified uses: | Metal-ceramic

1.4 Details of the manufacturer and importer:-

Manufacturer:

Ivoclar Vivadent Inc.
 175 Pineview Drive, Amherst, N.Y. 14228
 USA
 Tel. +1 800 533 6825
 Fax +1 716 691 2285

Ivoclar Vivadent Inc.
 1-6600 Dixie Road
 Mississauga, Ontario
 L5T 2Y2
 Canada
 Phone: +1 905 670 8499
 Fax: +1 905 670 3102

Importer:

Ivoclar Vivadent Ltd
 12 Omega St, Rosedale, Auckland, New Zealand
 Phone +64 9 914 9999 Fax +64 9 914 9990
 www.ivoclarvivadent.co.nz
 0800 764 766 (National Poison Centre)
 Poisons Hotline (24 hours / 7 days)

1.5 Emergency phone number:

2. HAZARD(S) IDENTIFICATION

2.1 GHS Classification:-

Void

2.2 GHS Label elements, including precautionary statements:-

Hazard Pictogram: | Void

Signal word: | Void

Hazard-determining components of labelling: | Void

Hazard statements: | Void

Precautionary statements: | Void

Additional information: | Dust generation
 Particular danger of slipping on leaked/spilled product.

2.3

Other hazards:-

Avoid breathing sanding dust /swarf.

Results of PBT and vPvB assessment;

PBT: | Not applicable.

vPvB: | Not applicable.

3 COMPOSITION/INFORMATION ON INGREDIENTS

Mixture of substances listed below with nonhazardous additions.

Ingredient name	CAS No.	Classification	Concentration
ceramic materials	66402-68-4	Void	Void

For the full text of the H-Statements mentioned in this Section, refer to Section 16.

4. FIRST AID MEASURES

4.1 Description of necessary first aid measures:-

If inhaled:

Supply fresh air; consult doctor in case of complaints.
In case of unconsciousness place patient stably in side position for transportation.

In case of skin contact:

Rinse with water.

In case of eye contact:

Generally the product does not irritate the skin.
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.

If swallowed:

Mechanical effects only.

Rinse out mouth and then drink plenty of water.

Symptoms caused by exposure:-

If symptoms persist consult doctor

4.2 Medical attention and special treatment:-

No further relevant information available.

4.3

No further relevant information available.

FIRE FIGHTING MEASURES

5. Suitable extinguishing equipment:-

5.1 Suitable extinguishing media:

Unsuitable extinguishing media:

Specific hazards arising from the substance/mixture/product:-

The product is not flammable.

The product is not flammable

5.2 Special protective equipment and precautions for fire fighters:-

No further relevant information available

5.3

Special personal protective equipment:

No further relevant information available

Precautions:

No further relevant information available

6. ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures:-

Particular danger of slipping on leaked/spilled product.

6.2 Environmental precautions:-

Do not allow to enter sewers/ surface or ground water.

6.3 Methods and materials for containment and cleaning up:-

Pick up mechanically.

Ensure adequate ventilation

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7. HANDLING AND STORAGE

7.1 Precautions for safe handling:-

Only adequately trained personnel should handle this product.

Prevent formation of dust.

Ensure good ventilation/exhaustion at the workplace.

7.2 Conditions for safe storage, including any incompatibilities:-

Store only in the original receptacle

Keep container tightly sealed.
Store in dry conditions.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Exposure control measures:-

Occupational exposure limits:

Component	CAS No.	Value	Parameters	Basis
Void	Void	Void	Void	Void

Ingredients with biological limit values:

It is recommended to observe the general dust limit of 10 mg/m³ (respirable fraction) / 3 mg/m³ (alveolar respirable fraction). The lists valid during the making were used as basis.

8.2 Biological monitoring:-

The lists valid during the making were used as basis.

Exposure controls / Personal protective equipment / General protective and hygienic measures:

Usual hygienic measures for dental practice.

Wash hands before breaks and at the end of work.

Keep away from foodstuffs, beverages and feed.

Vacuum clean contaminated clothing. Do not blow or brush off contamination.

Do not inhale dust / smoke / mist.

8.3 Control banding:-

Use good industrial hygiene practice and general ventilation.

8.4 Engineering controls:-

Void

8.5 Individual protection measures include PPE:-

Eye/face protection:



Safety glasses

Required

Skin protection:



Not required.

Respiratory protection:

Breathing equipment:

In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.

Recommended filter device for short term use: Filter P2

9. PHYSICAL/CHEMICAL PROPERTIES

9.1 Information on physical/chemical properties:-

a) **Appearance/Form:**

Powder

b) **Colour:**

According to product specification

c) **Odour:**

Odourless

d) **Odour threshold:**

Not determined.

e) **pH value:**

Not applicable.

f)	Melting point/melting range:	Undetermined.
g)	Boiling point/boiling range:	Not determined.
h)	Flash point:	Undetermined.
i)	Ignition temperature:	Undetermined.
j)	Self-igniting:	Product is not self-igniting.
k)	Danger of explosion:	Not determined.
l)	Upper/lower flammability or explosive limits:	Lower Not determined.
		Upper Not determined.
m)	Vapour pressure 20°C:	Not determined.
n)	Density at 20°C:	Not applicable.
o)	Relative density:	Not determined.
p)	Vapour density:	Not applicable.
q)	Evaporation rate:	Not applicable.
r)	Solubility in/miscibility with water:	Insoluble.
s)	Partition coefficient: n- octanol/water:	Not determined.
t)	Viscosity:	Dynamic Not applicable.
		Kinematic Not applicable.
u)	Solids content:	Not applicable.

10. STABILITY AND REACTIVITY

10.1 Reactivity:-

No further relevant information available.

10.2 Chemical stability:-

Stable under normal handling and storage conditions.

Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

10.3 Possibility of hazardous reactions:-

No dangerous reactions known.

10.4 Conditions to avoid:-

No further relevant information available.

10.5 Incompatible materials:-

No further relevant information available.

10.6 Hazardous decomposition products:-

None under normal conditions of storage and use.

11. TOXICOLOGICAL INFORMATION

11.1 Information on toxicological effects:-

Acute toxicity / Values relevant for classification:

LD/LC50 values relevant for classification:

107-88-0 butane-1,3-diol

Oral LD50 18610 mg/kg (rat)

Skin corrosion/irritation:

No irritant effect.

Serious eye damage/eye irritation:

Mechanical effects only.

Respiratory or skin sensitization:

No sensitising effects known.

Germ cell mutagenicity:

No further relevant information available.

Carcinogenicity:

No further relevant information available.

Reproductive toxicity:

No further relevant information available.

Specific target organ toxicity - single exposure:

No further relevant information available.

Specific target organ toxicity - repeated exposure:

No further relevant information available.

Aspiration hazard:

No further relevant information available.

Additional information:

No further relevant information available.

11.2 Information on possible routes of exposure:-

Short Term (Acute) Exposure:

No further relevant information available.

<p>Swallowed: Eyes: Skin: Inhaled: Long Term (Chronic) Exposure: Swallowed: Eyes: Skin: Inhaled:</p>	<p>No further relevant information available. No further relevant information available. No further relevant information available. No further relevant information available. No further relevant information available. No further relevant information available. No further relevant information available. No further relevant information available. No further relevant information available.</p>
<p>11.3 Early onset symptoms related to exposure:-</p>	<p>No further relevant information available.</p>
<p>11.4 Delayed health effects from exposure:-</p>	<p>No further relevant information available.</p>
<p>11.5 Exposure levels and health effects:-</p>	<p>No further relevant information available.</p>
<p>11.6 Interactive effects:-</p>	<p>No further relevant information available.</p>
<p>11.7 Other:-</p>	<p>No further relevant information available.</p>
<p>12. ECOLOGICAL INFORMATION</p>	
<p>12.1 Ecotoxicity:-</p>	<p>No further relevant information available.</p>
<p>12.2 Persistence/degradability:-</p>	<p>No further relevant information available.</p>
<p>12.3 Bioaccumulative potential:-</p>	<p>No further relevant information available.</p>
<p>12.4 Mobility in soil:-</p>	<p>No further relevant information available.</p>
<p>12.5 Other adverse effects:-</p>	<p>Generally not hazardous for water</p>
<p>Additional ecological information / General notes:-</p>	<p>No further relevant information available.</p>
<p>12.6 Other adverse effects:-</p>	<p>No further relevant information available.</p>
<p>13. DISPOSAL CONSIDERATIONS</p>	
<p>13.1 Disposal methods:-</p>	<p>Take to an approved landfill or a waste incineration plant, under conditions approved by the local authority. Disposal must be made according to official regulations.</p>
<p>14. TRANSPORT INFORMATION</p>	
<p>UN number ADR / IMDG / IATA:-</p>	<p>Void.</p>
<p>UN proper shipping name or technical name:-</p>	<p>Void.</p>
<p>ADR:</p>	<p>Void.</p>
<p>IMDG, IATA:</p>	<p>Void.</p>
<p>Transport hazard class(es):</p>	<p>Void.</p>
<p>Label:</p>	<p>Void.</p>
<p>Packaging group:</p>	<p>Void.</p>
<p>Environmental hazards:</p>	<p>Not applicable.</p>
<p>Special precautions for user:</p>	<p>Not applicable.</p>
<p>Danger code:</p>	<p>Void.</p>
<p>EMS Number:</p>	<p>Void.</p>
<p>Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:</p>	<p>Not applicable.</p>
<p>Transport/Additional information:</p>	<p>Product is not classified as a dangerous good for transport (ADR, IMDG, IATA).</p>

15. REGULATORY INFORMATION

15.1 Safety, health and environmental regulations/legislation specific for the substance/mixture/product:-

Sara

- Section 355 (extremely hazardous substances):
Substance is not listed.
- Section 313 (Specific toxic chemical listings):
Substance is not listed.
- TSCA (Toxic Substances Control Act):
Substance is listed.

· **Proposition 65**

- Chemicals known to cause cancer:
Substance is not listed.
- Chemicals known to cause reproductive toxicity for females:
Substance is not listed.
- Chemicals known to cause reproductive toxicity for males:
Substance is not listed.
- Chemicals known to cause developmental toxicity:
Substance is not listed.

Carcinogenic categories

- EPA (Environmental Protection Agency)
Substance is not listed.
- TLV (Threshold Limit Value established by ACGIH)
Substance is not listed.
- NIOSH-Ca (National Institute for Occupational Safety and Health)
Substance is not listed.
- OSHA-Ca (Occupational Safety & Health Administration)
Substance is not listed.

15.2 Chemical safety assessment:

A Chemical Safety Assessment has not been carried out.
This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

16. OTHER INFORMATION

Key to abbreviations/acronyms used in SDS:- | Void

Key literature references/data sources used to compile SDS:-

Standard EN420:2003 General requirements for protective gloves: disposable gloves, e.g. nitrile rubber, material thickness 0.1 mm (Australian Standard 2161).
Long-term exposure (Level 6: < 480 min): protective gloves, e.g. nitrile rubber, material thickness 0.7 mm (Australian Standard 2161).
Personal eye protection - Eye and face protectors for occupational applications: safety glasses (Australian Standard AS 1336 and AS/NZS 1337.1:2010).

Copyright statement:-

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Abbreviations and acronyms:-

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road).

IMDG: International Maritime Code for Dangerous Goods.

IATA: International Air Transport Association.

GHS: Globally Harmonised System of Classification and Labelling of Chemicals.

EINECS: European Inventory of Existing Commercial Chemical Substances.

ELINCS: European List of Notified Chemical Substances.

CAS: Chemical Abstracts Service (division of the American Chemical Society).

LC50: Lethal concentration, 50 percent.

LD50: Lethal dose, 50 percent.

Flam. Liq. 2: Flammable liquids, Hazard Category 2.

Flam. Liq. 3: Flammable liquids, Hazard Category 3.

Acute Tox. 4: Acute toxicity, Hazard Category 4.

Skin Irrit. 2: Skin corrosion/irritation, Hazard Category 2.

Eye Irrit. 2: Serious eye damage/eye irritation, Hazard Category 2.

Repr. 2: Reproductive toxicity, Hazard Category 2.

STOT SE 3: Specific target organ toxicity - Single exposure, Hazard Category 3.

STOT RE 2: Specific target organ toxicity - Repeated exposure, Hazard Category 2.

Asp. Tox. 1: Aspiration hazard, Hazard Category 1.

Disclaimer:-

Ivoclar Vivadent Ltd believes the information in this document to be accurate as at the date of preparation, but, to the maximum extent permitted by law, Ivoclar Vivadent Ltd accepts no responsibility for any loss or damage caused by any person acting or refraining from action because of this information.

The provision of this information should not be construed by anyone as a recommendation to use this product. In particular, no one should use any product in violation of any patent or other intellectual proprietary rights or in breach of any statute or regulation.

Users should rely on their own knowledge and inquiries and make their own determination as to the applicability of this information in relation to their particular purposes and specific circumstances. Each user should read this SDS and consider the information in the context of how the product will be handled and used in the workplace and in conjunction with other substances or products.

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.