

Programat® CS6



Sintering Furnace Accessories

The ideal accessories for your Programat® CS6

Sintering Furnace Temperature Checking Set Starter Kit

627380

The Temperature Checking Set Starter Kit comprises all the components necessary to check and calibrate the sintering temperature of the Programat CS6.



Programat® Calibration Measuring Aid

752035

The Programat Calibration Measuring Aid supports users in measuring the control rings and provides them with a reliable measuring result as a basis for the furnace calibration.



Sintering Furnace Temperature Checking Set UTH (880 °C) (Refill Package)

688632

10 sintering rings to check the exact crystallization temperature of the Programat CS6.



Sintering Furnace Temperature Checking Set MTH (1500 °C) (Refill Package)

627398

10 sintering rings to check the exact sintering temperature of the Programat CS6.



Programat® WLAN-Kit

655949

The WLAN Kit is used for the wireless data transfer.



Symbolic figure

Programat® CS6



Sintering Furnace Accessories

The ideal accessories for your Programat® CS6

Programat® CS6 Tweezer Holder

747590

The tweezer holder can be attached to the housing of the Programat CS6. It enables you to conveniently remove and replace the tweezers.



Programat® CS6 Tweezers

749897

The special tweezers with the ceramic tips are used to load the furnace.



Programat® CS6 Object Plate

746379

On the object plate of the Programat CS6, restorations are placed in the firing chamber. The plate is also used for cooling down zirconia restorations outside the firing chamber.



Digital Shade Assistant Starter Kit

679917

The Digital Shade Assistant Starter Kit is especially designed for digital shade selection. It includes all the components the user needs. The Digital Shade Assistant Holder is used for holding the shade guide teeth.


Delivery Form:
3 Digital Shade Assistant Holder
1 A-D shade guide



IPS Speed Tray

746355

The IPS Speed Tray is required for the crystallization and glazing programs of the Programat CS6.



Digital Shade Assistant Holder Kit

679916

The Digital Shade Assistant Holder is used for holding the shade guide teeth for the DSA function.

Delivery Form:
2 Digital Shade Assistant Holders

