



Skyce[®]
Tooth jewels

**A touch
of glamour**

ivoclar
vivadent[®]
passion vision innovation

Create sparkling accents

The dazzling Skyce® crystals add a flash of extra sparkle to a smile. The tooth jewels come with a coordinated adhesive system so that you can provide the high-quality professional results that your patients expect.

Benefits for your customers:

Preservation of tooth structure due to non-invasive preparation and bonding procedure

Coordinated product system with Total Etch and Heliobond to ensure a reliable bond and a gentle residue-free removal

Available in two sizes:
 ø 1.9 mm and ø 2.4 mm

Magical light effects
 due to premium quality crystal glass



Etch dental enamel with phosphoric acid



Apply Heliobond adhesive



Place Skyce® onto the bonding site



Light cure from all sides

Source:
 Dr A. Peschke, Ivoclar Vivadent AG,
 Schaan, Liechtenstein

Art. no.	Description
559320AN	Skyce® Intro Package
559322AN	Skyce® Refill 5 Crystal, large DE
559321AN	Skyce® Refill 5 Crystal, small DE



Production and Distribution
Ivoclar Vivadent AG
 Benderstr. 2
 9494 Schaan
 Liechtenstein
 Tel. +423 235 35 35
 Fax +423 235 33 60
www.ivoclarvivadent.com

ivoclar
vivadent
 passion vision innovation