

## Combination table

IPS Style Ceram																									
A-D	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4					
IPS Style Ceram Opaquer	0 BL1/BL2		0 BL3/BL4		0 A1	0 A2	0 A3	0 A3.5	0 A4	0 B1	0 B2	0 B3	0 B4	0 C1	0 C2	0 C3	0 C4	0 D2	0 D3	0 D4					
IPS Style Ceram Intensive Opaquer					IO white				IO violet	IO brown				IO incisal											
IPS Style Ceram Margin	M BL				M 1	M 2	M 3	M 4	M 7	M 1	M 2	M 4	M 4	M 5	M 6	M 6	M 7	M 5	M 3	M 6					
IPS Style Ceram Deep Dentin	DD BL1/BL2		DD BL3/BL4		DD A1	DD A2	DD A3	DD A3.5	DD A4	DD B1	DD B2	DD B3	DD B4	DD C1	DD C2	DD C3	DD C4	DD D2	DD D3	DD D4					
IPS Style Ceram Dentin	D BL1	D BL2	D BL3	D BL4	D A1	D A2	D A3	D A3.5	D A4	D B1	D B2	D B3	D B4	D C1	D C2	D C3	D C4	D D2	D D3	D D4					
IPS Style Ceram Incisal	I BL				I 1	I 2	I 3	I 4	I 5	I 1	I 2	I 3	I 4	I 2	I 3	I 4	I 5	I 2	I 3	I 4					
IPS Style Ceram Add-On	A-0 BL				A-0 690°				A-0 Incisal		A-0 Dentin		A-0 Margin												
IPS Style Ceram Impulse	Occlusal Dentin	OD orange	OD brown	Mamelon				MM light	MM yellow-orange	MM salmon	Opal Effect				OE 1	OE 2	OE 3	OE 4	OE 5	OE violet	Incisal Edge	Incisal Edge			
	Transpa	T neutral	T clear	T blue	T brown-grey	T orange-grey	Special Incisal		SI yellow	SI grey	Inter Incisal		II white-blue		Cervical Transpa				CT yellow	CT orange-pink	CT khaki	CT orange			
IPS Style Ceram Gingiva	Gingiva Opaquer	G0 pink	Basic Gingiva				BG34	Gingiva				G1	G2	G3	G4	G5	Intensive Gingiva				IG1	IG2	IG3	IG4	IG5
IPS Ivocolor Essence	E01 white	E02 creme	E03 lemon	E04 sunset	E05 copper	E06 hazel	E07 olive	E08 khaki	E09 terracotta	E10 mahogany	E11 cappuccino	E12 espresso	E13 terra	E14 profundo	E15 ocean	E16 sapphire	E17 anthracite	E18 black	E19 rose	E20 coral					
									E21 basic red		E22 basic yellow		E23 basic blue												
IPS Ivocolor Shade	Shade 0				Shade 1		Shade 2		Shade 3		Shade 4		Shade 5		Shade 6				Shade 7		Shade 6				
									Shade Incisal 1		Shade Incisal 2		Shade Incisal 3												
IPS Style Ceram One																									
A-D	BL1	BL2	BL3	BL4	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4					
IPS Style Ceram Opaquer	IO white	0 BL1/BL2	0 BL3/BL4		0 A1	0 A2	0 A3	0 A3.5	0 A4	0 B1	0 B2	0 B3	0 B4	0 C1	0 C2	0 C3	0 C4	0 D2	0 D3	0 D4					
IPS Style Ceram One	One BL				One 1	One 2	One 3	One 4	One 7	One 1	One 2	One 4	One 4	One 5	One 6	One 6	One 7	One 5	One 3	One 6					